

June 2020

Dear Assessor,

Please read the following information regarding changes to the ACC MRES Simple List Equipment Walking Frames which take effect Wednesday 1st July 2020.

There is an ISO standard being introduced for Walking Frames that will impact the Simple List Equipment. This ISO standard applies for all new equipment supplied going forward only and does not apply retrospectively to equipment already with clients.

- ISO 11199-2:2005 Walking aids manipulated by both arms -- Requirements and test methods -- Part 2: Rollators

Any Walking Frames with wheels of less than 180mm (7 inches) in diameter, are now considered suitable for indoor use **only**.

As a result, we have had to respond to this requirement by:

- introducing 2 new Walking Frames that have 8" Wheels that **are** suitable for **both** indoor and outdoor use
- remove 2 current Walking Frames that have 7" Wheels that **are not** suitable for **both** indoor and outdoor use
- inform you that 8 current Walking Frames are now **only** suitable for indoor use

Please take note of the following SIC code and description updates when requesting and ordering the newly introduced options suitable for indoor and outdoor use.

**New Simple List Equipment Walking Frames with 8" Wheels**

Simple List SIC Code	Equipment Description
EM036N	Walking Frame 8" Wheels <b>with</b> Tray indoor & outdoor use
EM035N	Walking Frame 8" Wheels <b>without</b> Tray indoor & outdoor use

## Removal of Simple List Equipment Walking Frames with 7” Wheels

Simple List SIC Code	Equipment Description
EM028N	Walking Frame, Super Stroller <b>with</b> Tray
EM038N	Walking Frame, Super Stroller <b>without</b> Tray

## Updated Walking Frames ONLY Suitable for Indoor Use

Simple List SIC Code	Equipment Description
EM034N	Pacer Front Wheels Rear Stops Basket and Tray ( <i>Only Suitable for Indoor Use</i> )
EM031N	Pacer Castors Rear Stops Basket and Tray ( <i>Only Suitable for Indoor Use</i> )
EM039N	Pacer Front Castors, Rear Wheels, Handbrakes, Basket & Tray ( <i>Only Suitable for Indoor Use</i> )
EM033N	Pacer Front Castors Rear Wheels, Basket and Tray ( <i>Only Suitable for Indoor Use</i> )
EM032N	Pacer Front Castors, Rear Skis, Basket & Tray ( <i>Only Suitable for Indoor Use</i> )
EM029N	Pacer Gutter Walker Front Castors, Rear Wheels, Handbrakes ( <i>Only Suitable for Indoor Use</i> )
EM030N	Pacer Gutter Front Castors, Rear Wheels ( <i>Only Suitable for Indoor Use</i> )
EM037N	Super Stroller Narrow Walking Frame ( <i>Only Suitable for Indoor Use</i> )

These changes take effect on Wednesday 1st July 2020 and retrospective replacement for Walking Frames currently with clients is not required and will not apply.

## Connecting with Us

For any enquires you may have, please contact us in the usual manner.

Website: [www.accessable.co.nz](http://www.accessable.co.nz)

Email: [acc@accessable.co.nz](mailto:acc@accessable.co.nz)

Freephone: 0508 001 002

Sapphire Lautua  
ACC Contract Relationship Manager