



# Accessible Power Wheelchairs Product Overview and Resource

**May 2025**

# Table of Contents



Content	Page
Accessible Power Wheelchairs Categorisation Framework	4
Power Wheelchair Products List	6
aSAP Catalogue – Searching Tips	7
aSAP Catalogue Product Details	12
Accessible List Paediatric Power Wheelchair Range	15
- Paediatric Power Wheelchair Products List Overview	16
- Paediatric Power Wheelchair with Moderate Configuration	18
- Paediatric Power Wheelchair with High Configuration	27
Accessible List Adult Power Wheelchair Range	34
- Adult Power Wheelchair Products List Overview	35
- Differences between power wheelchair front, mid and rear drive wheel type	38

# Table of Contents



Content	Page
- Adult Power Wheelchair with Minimal Configuration	39
- Adult Power Wheelchair with Moderate Configuration	46
- Adult Power Wheelchair with High Configuration	67
Accessible List Bariatric Power Wheelchair Range	86
- Bariatric Power Wheelchair Products List Overview	87
- Bariatric Power Wheelchair with Moderate Configuration	89
- Bariatric Power Wheelchair with High Configuration	96
Acknowledgement	103

# Accessible Power Wheelchairs Categorisation Framework



Accessible's power wheelchair categorisation framework is designed to deliver a comprehensive selection of power wheelchairs that meet the diverse needs of our clients. We provide options ranging from entry-level wheelchairs to those suitable for clients with high and complex needs. Our offerings also include power wheelchairs specifically designed for paediatric, adult, and bariatric clients, ensuring that all users can find a suitable option for their mobility requirements.

Power Wheelchair Categories	Power Wheelchair with Minimal Configuration	Power Wheelchair with Moderate Configuration	Power Wheelchair with High Configuration
Paediatric	User Weight: Up to 136kg Seat Width: Up to 16"		
Adult	User Weight: Up to 175kg Seat Width: Up to 22"		
Bariatric	User Weight: Greater than 175kg Seat Width: Greater than 22"		



Power Wheelchair Categories	Power Wheelchair with Minimal Configuration	Power Wheelchair with Moderate Configuration	Power Wheelchair with High Configuration
<b>Clinical Profile Description</b>			
<b>Profile</b>	<ul style="list-style-type: none"> <li>• Suitable for indoor mobility and immediate outdoor environment around home</li> <li>• Individuals new to using power wheelchairs</li> <li>• Can support individuals displaying WMPM (Wheeled Mobility and Postural Management) Level 1 indicators</li> </ul>	<ul style="list-style-type: none"> <li>• Suitable for indoor mobility and outdoor use in urban environments</li> <li>• Designed for individuals with moderate to complex needs, presenting a stable condition</li> <li>• Appropriate for a majority of clients with various conditions</li> <li>• Can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1)</li> </ul>	<ul style="list-style-type: none"> <li>• Suitable for indoor mobility and outdoor urban terrain</li> <li>• Designed for individuals with high and complex needs, including progressive conditions with the potential for changing needs</li> <li>• Features multiple power seat functions</li> <li>• Supports alternative controls</li> <li>• Can support individuals displaying WMPM Level 2 clinical indicators</li> </ul>
<b>Standard Features</b>	<ul style="list-style-type: none"> <li>• Controller (non-expandable electronics)</li> <li>• No power seat functions</li> <li>• Standard arm/ leg/ foot supports</li> <li>• Standard hip belt</li> <li>• Battery charger</li> </ul>	<ul style="list-style-type: none"> <li>• Controller (up to 2 functions/ expandable electronics)</li> <li>• Maximum of 2 power seat functions, including power tilt</li> <li>• Adjustable seat pan</li> <li>• Adjustable back cane angle</li> <li>• Arm/ leg/ foot supports</li> <li>• Off-the-shelf seating compatible</li> <li>• Hip belt</li> <li>• Battery charger</li> <li>• Transit tie down</li> </ul>	<ul style="list-style-type: none"> <li>• Controller (expandable electronics)</li> <li>• More than 2 power seat functions</li> <li>• Adjustable seat pan</li> <li>• Adjustable back cane angle</li> <li>• Arm/ leg/ foot supports</li> <li>• Off-the-shelf seating compatible</li> <li>• Hip belt</li> <li>• Battery charger</li> <li>• Transit tie down</li> </ul>

# Power Wheelchair Products List



Power Wheelchair Categories	Power Wheelchair with Minimal Configuration	Power Wheelchair with Moderate Configuration	Power Wheelchair with High Configuration
Paediatric	-	<b>Z500 Mini (Medifab)</b> Quantum Edge 3 Stretto (Allied Medical) Alltrack P3 (Cubro) TA IQ Junior (Cubro)	Quantum Edge 3 Stretto (Allied Medical) Alltrack P3 (Cubro) TA IQ Junior (Cubro)
Adult	<u>MWD options:</u> <b>M1 (non-tilt) (Permobil)</b>  <u>RWD options:</u> <b>Karma KP25 2.0 (Allied Medical)</b> <b>Q100 RWD (Medifab)</b> <b>Q200 RWD (Medifab)</b>	<u>MWD options:</u> <b>M1 (with tilt) (Permobil)</b> <b>Pride J623 2.0 (Allied Medical)</b> <b>Q500M (Medifab)</b> TA IQ MWD (Cubro) TDXSP2 REM216 (Invacare) Quantum Edge 3 (Allied Medical) M3 Corpus (Permobil)  Q500 Mini (Medifab) Quantum Edge 3 Stretto (Allied Medical)  <u>RWD options:</u> Aviva RX REM216 (Invacare)	<u>MWD options:</u> <b>Q500M (Medifab)</b> TA IQ MWD (Cubro) TDXSP2 REM400 (Invacare) Quantum Edge 3 (Allied Medical) M3 Corpus (Permobil) Q500 Mini (Medifab) Quantum Edge 3 Stretto (Allied Medical)  <u>RWD options:</u> Aviva RX REM400 (Invacare)  <u>FWD options:</u> <b>F3 Corpus (Permobil)</b>
Bariatric	-	<u>MWD options: (Up to 204kg SWL)</u> TDXSP2 HD REM216 (Invacare) Quantum Q6 Edge HD (Allied Medical)  <u>FWD option: (&gt;204kg, up to 272kg)</u> Quantum Q1450 (Allied Medical)	<u>MWD options: (Up to 204kg SWL)</u> TDXSP2 HD REM400 (Invacare) Quantum Q6 Edge HD (Allied Medical)  <u>FWD option: (&gt;204kg, up to 272kg)</u> Quantum Q1450 (Allied Medical)

**Green** = power wheelchairs in one category only

**Blue** = power wheelchairs in more than one category

MWD = Mid Wheel Drive

RWD = Rear Wheel Drive

FWD = Front Wheel Drive

# aSAP Catalogue – Searching Tips

# aSAP Catalogue Tips for Searching

- Search by classification category

## Accessible aSAP Portal

Customers ▾ Advice Requests ▾ Applications ▾ Catalogue ▾ Accessible MOH Assessor ▾

asap Portal > Full Catalogue (MOH)

### Full Catalogue (MOH)

Access Equipment

Bathroom and Toilet

Bedroom

Communication Assistive Technology

Hearing Assistive Technology (HAT)

Kitchen, Dining, Living, Household Management

**Mobility**

Postural Management (Seating, Standing and Lying)

Removable Housing Equipment

Transfers and Manual Handling

BAND SELECTION

☐ Refurb

☐ Trial

SIC code

Supplier Code

Search By Supplier

Product Name Search:

Search By Keyword 1

Search By Keyword 2

Tips for searching the catalogue:

1. Select the product group by using the blue buttons on the left-hand side.

2. Refine the search with the filters above. If using keywords for product searches, these need to be entered in the search box.

3. Then push

Wheelchair - Power Front Wheel Drive

Wheelchair - Power Front Wheel Drive Multi

**Wheelchair - Power Mid Wheel Drive**

Wheelchair - Power Mid Wheel Drive Multi

Wheelchair - Power Rear Wheel Drive

Wheelchair - Power Rear Wheel Drive Multi

Item Name	SIC	Supplier Code	Price	Supplier	Band Type
Alltrack P3 12-16inWx12-18inD Tilt Elevate Powerchair	EM323	0600-0026		Cubro	Band 3
Alltrack P3 12-16inWx12-18inD Tilt Powerchair	EM322	0600-0024		Cubro	Band 3
Alltrack P3 12-16inWx12-18inD Tilt Recline Powerchair	EM324	0600-0025		Cubro	Band 3
Alltrack P3 12-16inWx12-18inD TRE (PCM) Powerchair	EM326	0600-0028		Cubro	Band 3
Alltrack P3 12-16inWx12-18inD TRE Powerchair	EM325	0600-0027		Cubro	Band 3
Edge 3 MWD 12-16inWx12-16inD Tilt Recline Powerchair	EM378	PKE3TRFC12161216		Allied Medical Ltd	Band 3
Edge 3 MWD 12-16inWx16-20inD Tilt Powerchair	EM376	PKE3TFC12161620		Allied Medical Ltd	Band 3
Edge 3 MWD 12-16inWx16-20inD Tilt Recline Powerchair	EM378	PKE3TRFC12161620		Allied Medical Ltd	Band 3
Edge 3 MWD 14-16inWx16-20inD Tilt Elevate Powerchair	EM377	PKE3TFC20242024		Allied Medical Ltd	Band 3
Edge 3 MWD 14-16inWx16-20inD TRE (PCM) Powerchair	EM381	PKE3TERFPC14161620		Allied Medical Ltd	Band 3
Edge 3 MWD 14-16inWx16-20inD TRE Powerchair	EM379	PKE3TERFPC14161620		Allied Medical Ltd	Band 3
Edge 3 MWD 16-20inWx12-16inD Tilt Powerchair	EM376	PKE3TFC16201216		Allied Medical Ltd	Band 3
Edge 3 MWD 16-20inWx12-16inD Tilt Recline (PCM) Powerchair	EM380	PKE3TRFPC16201216		Allied Medical Ltd	Band 3
Edge 3 MWD 16-20inWx12-16inD Tilt Recline Powerchair	EM378	PKE3TRFPC16201216		Allied Medical Ltd	Band 3
Edge 3 MWD 16-20inWx16-20inD Tilt Recline (PCM) Powerchair	EM380	PKE3TRFPC16201620		Allied Medical Ltd	Band 3
Edge 3 MWD 16-20inWx16-20inD TRE (PCM) Powerchair	EM381	PKE3TERFPC16201620		Allied Medical Ltd	Band 3
Edge 3 MWD 16-20inWx16-20inD TRE Powerchair	EM379	PKE3TERFPC16201620		Allied Medical Ltd	Band 3



# aSAP Catalogue Tips for Searching

- All power wheelchairs are Band 3 items
  - **Apply Band 3 filter** to concentrate search for power wheelchairs

- Search by **Product Name**

The screenshot displays the aSAP Catalogue search interface. A 'BAND SELECTION' dropdown menu is open, showing four options: 'Band 1 List', 'Band 2 List', 'Band 3 List' (which is selected and highlighted with a blue border), and 'Non List'. Below the dropdown is a green 'Apply' button. In the background, the main search interface is visible, featuring a 'BAND SELECTION' button, a 'Refurb' checkbox, a 'Trial' checkbox, a 'SIC code' input field, a 'Supplier Code' input field, and a 'Product Name Search:' section with two input fields labeled 'Search By Keyword 1' and 'Search By Keyword 2'. There are also buttons for 'Search By Supplier', 'Search By Keyword 1', and 'Search By Keyword 2'. At the bottom right, there are 'SEARCH' and 'RESET SEARCH' buttons.

- Use both fields
- Use the keyword “**Powerchair**”, not “Power wheelchairs”

# aSAP Catalogue Tips for Searching

- **Product Name Search** – search by power seat functions required

- Power tilt – “tilt powerchair”

Product Name Search:

- Power tilt and recline

Product Name Search:

- Power tilt and elevate

Product Name Search:

- Power tilt, recline and seat elevate (TRE)

Product Name Search:

- Power tilt and powered centre-mount (PCM) elevating leg rest

- Centre-mount (CM)


Product Name Search:

- Power tilt and powered swing away (PSA) elevating leg rest

- Swing away (SA)

Product Name Search:

# aSAP Catalogue Tips for Searching

- **Note:** the majority of power wheelchair packages come with a centre-mount leg rest configuration, and some come standard with a swing-away style. Please click on the  icon of the product and refer to the **Key feature** section.

- Example:

Item Name	SIC	Supplier Code	Price	Supplier	Band Type	Item Type
Alltrack P3 12-16inWx12-18inD Tilt Elevate Powerchair	EM323	0600-0026		Cubro	Band 3	List



## Features

## Specifications

## Support

## Accreditations

## Classifications

### Intended Use

Power wheelchair with moderate configuration (Paediatric, Adult). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 2 clinical indicators.

### Description

Alltrack P3 MWD power wheelchair - Tilt/elevate. Base width: 535mm. Base length: 890mm. Seat to floor height: 445mm. Turning radius: 510mm. Speed: 10km/h. Battery capacity: 50amps. Tilt 50°. Seat Elevate 11in.

### Key Features

Package includes: Standard controller and electronics (R-net expandable). Swing away controller mount / receiver. Centre mount leg hanger with flip up, angle adjustable footplates (2 individual pieces). Flip back, height adjustable arm supports with standard type, full length arm pads. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery charger.

- **Viewing details – Example view**



## 12

# aSAP Catalogue Product Details

## Intended Use:

- Provides information on the **power wheelchair category**
- **Intended type of users** and **environment** of use

Level of Configuration
Minimal configuration
Moderate configuration
High configuration

Type of Users
Paediatric
Adult
Bariatric

- Considerations with reference to the **Wheeled Mobility and Postural Management (WMPM) indicators** – Level 1 and/or 2
- Example: **Intended Use**

Power wheelchair with moderate configuration (Paediatric, Adult). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 2 clinical indicators.

# aSAP Catalogue Product Details

## Description:

- Provides the **Product Name** and **Key Functions** (including **powered functions**) of the product.
- Example: **Description**

Alltrack P3 MWD power wheelchair - Tilt/elevate. Base width: 535mm. Base length: 890mm. Seat to floor height: 445mm. Turning radius: 510mm. Speed: 10km/h. Battery capacity: 50amps. Tilt 50°. Seat Elevate 11in.

## Key Features:

- Provides more detailed information about **items included in the package** (e.g., controller, arm and leg supports, accessories)
- Example: **Key Features**

Package includes: Standard controller and electronics (R-net expandable). Swing away controller mount / receiver. Centre mount leg hanger with flip up, angle adjustable footplates (2 individual pieces). Flip back, height adjustable arm supports with standard type, full length arm pads. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery charger.

# Accessible List Paediatric Power Wheelchair Range



# Paediatric Power Wheelchair Products List Overview



Power Wheelchair Categories	Power Wheelchair with Minimal Configuration	Power Wheelchair with Moderate Configuration	Power Wheelchair with High Configuration
Paediatric	-	<p><b>Z500 Mini (Medifab)</b></p> <p><b>Quantum Edge 3 Stretto (Allied Medical)</b></p> <p><b>Alltrack P3 (Cubro)</b></p> <p><b>TA IQ Junior (Cubro)</b></p>	<p><b>Quantum Edge 3 Stretto (Allied Medical)</b></p> <p><b>Alltrack P3 (Cubro)</b></p> <p><b>TA IQ Junior (Cubro)</b></p>

**Green** = power wheelchairs in one category only

**Blue** = power wheelchairs in more than one category

MWD = Mid Wheel Drive

RWD = Rear Wheel Drive

FWD = Front Wheel Drive



# Paediatric Power Wheelchair Range

Within this range, there are power wheelchair options with:

- **Seat dimensions with growth availability** that cater to the sizing of very young clients, as well as those in their teenage and young adult years
- **Minimal configuration through to highly configurable options** that meet clients' needs across the range of complexity
- A **compact base and a small turning radius** for use in indoor environments with confined spaces, e.g. school/ classroom environment
- **Seat functions**, ranging from tilt options to multiple power seat functions
- **Controller capabilities** with expandable electronics
- **Varying degrees of motor/ battery capacities**

*Note: When selecting power wheelchairs and features, it's important to consider funding criteria. Assessors should determine the essential needs based on the client's context, such as the necessity for a power wheelchair with a compact base, considering the client's environment and the requirements for essential access.*

# Paediatric Power Wheelchair with Moderate Configuration



Power Wheelchair Categories	Power Wheelchair with Moderate Configuration
Paediatric	<p><b>Z500 Mini (Medifab)</b></p> <p><b>Quantum Edge 3 Stretto (Allied Medical)</b></p> <p><b>Alltrack P3 (Cubro)</b></p> <p><b>TA IQ Junior (Cubro)</b></p>

**Green** = power wheelchairs in one category only

**Blue** = power wheelchairs in more than one category

MWD = Mid Wheel Drive

RWD = Rear Wheel Drive

FWD = Front Wheel Drive

# Z500 Mini MWD Powerchair – Moderate Configuration



**Description:** Z500 Mini MWD Paediatric Powerchair, Power Tilt, Elevate, Adjustable Seat, Armrests & Footrests. 13in Drive Wheels. 7in Castors. 50 Ah Batteries. 10km Speed. 30km Range. 50 Degree Posterior Tilt. 30cm Seat Elevate. 40 Degree Manual Recline. Manual swing away (SA) Removable Legrests, swings In/Out. Angle Adjustable Footplates.

**Cost bracket:** \$16000.00 - \$17000 (Tilt Elevate)

**Supplier:** Medifab Limited

**Intended Use:** Power wheelchair with moderate configuration (Paediatric). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators.

**Key Features:** Package includes: R-Net Expandable Electronics. 13in Drive Wheels. 7in Castors. Transit Tie Downs. 4 Pole Motors. 50Ah Gel Batteries. Battery Charger. ASAP Adjustable 10-14in or 12-16in Seat Top. Power Seat Tilt. Power Seat Elevate. Manual Recline. 80 Degree Swing Away Legrests. Aluminium Angle Adjustable Footplates. Cantilever Flip Back Armrests. RHS Swing Away Controller Mount. CJSM1 Standard Joystick. Standard Lap Belt. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg

# Z500 Mini MWD Powerchair – Moderate Configuration



**Description:** Z500 Mini MWD Paediatric Powerchair, Power Tilt, Elevate, Adjustable Seat, Armrests & Footrests. 13in Drive Wheels. 7in Castors. 50 Ah Batteries. 10km Speed. 30km Range. 50 Degree Posterior Tilt. 30cm Seat Elevate. 40 Degree Manual Recline. Manual swing away (SA) Removable Legrests, swings In/Out. Angle Adjustable Footplates.

**Cost bracket:** **\$16000.00 - \$17000 (Tilt Elevate)**

**Supplier:** Medifab Limited

## Product Items on aSAP

Item Name	SIC	Supplier Code
Z500 Mini MWD 10-14inWx10-16inD Tilt Elevate (SA) Powerchair	EM363	1660-1901-001
Z500 Mini MWD 12-16inWx12-18inD Tilt Elevate (SA) Powerchair	EM364	1660-1901-001

## Product URL

<https://medifab.com/nz/product/zippie-z500m-mini/>

# Stretto MWD Powerchair – Moderate Configuration



**Description:** The narrowest, most compact and manoeuvrable power base in the Quantum range. 14in Drive Wheels. 5in Solid Castor Wheels Front. 6in Solid Castor Wheels Rear. Min Seat to Floor Height 17.5in. Maximum Rate Slope 6°. Max Speed 10km/h. Mid-Wheel Drive (MWD). Q-Logic 3 Expandable Electronics. SRS Suspension. High Efficiency And Mileage. 2 Pole Motor. Fender lights. Selections of: Power tilt. Power tilt and elevate. Power Tilt and Recline.

**Cost bracket:** \$15000.00 - \$16000 (Tilt), \$18000.00 - \$19000 (Tilt Elevate), \$19000.00 - \$20000 (Tilt Recline)

**Supplier:** Allied Medical Ltd

**Intended Use:** Power wheelchair with moderate configuration (Paediatric, Adult). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1).

**Key Features:** Package includes: Standard controller and electronic (Q-Logic 3 expandable). Swing away controller mount / receiver. Centre mount leg hanger with flip up, angle adjustable footplates. Flip back, height adjustable arm supports. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery and charger. Selections of: Power tilt. Power tilt and elevate. Power Tilt and Recline. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg

# Stretto MWD Powerchair – Moderate Configuration



## Product Items on aSAP

Item Name	SIC	Supplier Code
Stretto MWD 12-16inWx12-16inD Tilt Powerchair	EM370	PKSTTFC12161216
Stretto MWD 12-16inWx16-20inD Tilt Powerchair	EM370	PKSTTFC12161620
Stretto MWD 16-20inWx12-16inD Tilt Powerchair	EM370	PKSTTFC16201216
Stretto MWD 16-20inWx16-20inD Tilt Powerchair	EM370	PKSTTFC16201620
Stretto MWD 16-20inWx20-22inD Tilt Powerchair	EM370	PKSTTFC16202024

## Product Items on aSAP

Item Name	SIC	Supplier Code
Stretto MWD 12-16inWx12-16inD Tilt Elevate Powerchair	EM371	PKSTTEFC12161216
Stretto MWD 14-16inWx16-20inD Tilt Elevate Powerchair	EM371	PKSTTEFC14161620
Stretto MWD 16-20inWx16-20inD Tilt Elevate Powerchair	EM371	PKSTTEFC16201620

Item Name	SIC	Supplier Code
Stretto MWD 12-16inWx12-16inD Tilt Recline Powerchair	EM372	PKSTTRFC12161216
Stretto MWD 12-16inWx16-20inD Tilt Recline Powerchair	EM372	PKSTTRFC12161620
Stretto MWD 16-20inWx12-16inD Tilt Recline Powerchair	EM372	PKSTTRFC16201216
Stretto MWD 16-20inWx16-20inD Tilt Recline Powerchair	EM372	PKSTTRFC16201620

## Product URL

<https://alliedmedical.co.nz/products/quantum-edge-3-stretto>



# Alltrack P3 MWD Powerchair – Moderate Configuration



**Description:** Alltrack P3 MWD power wheelchair. Base width: 535mm. Base length: 890mm. Seat to floor height: 445mm. Turning radius: 510mm. Speed: 10km/h. Battery capacity: 50amps. Selections of: Power Tilt 50°, Tilt and Seat elevate 11in, Tilt and recline with mechanical shear reduction seating functions.

**Cost bracket:** \$12000.00 - \$13000 (Tilt), \$15000.00 - \$16000 (Tilt Elevate), \$19000 - \$20000 (Tilt Recline)

**Supplier:** Cubro

**Intended Use:** Power wheelchair with moderate configuration (Paediatric). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1).

**Key Features:** Package includes: Standard controller and electronics (R-net expandable). Swing away controller mount / receiver. Centre mount leg hanger with flip up, angle adjustable footplates (2 individual pieces). Flip back, height adjustable arm supports with standard type, full length arm pads. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery charger. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg

# Alltrack P3 MWD Powerchair – Moderate Configuration



**Description:** Alltrack P3 MWD power wheelchair. Base width: 535mm. Base length: 890mm. Seat to floor height: 445mm. Turning radius: 510mm. Speed: 10km/h. Battery capacity: 50amps. Selections of: Power Tilt 50°, Tilt and Seat elevate 11in, Tilt and recline with mechanical shear reduction seating functions.

**Cost bracket:** \$12000.00 - \$13000 (Tilt), \$15000.00 - \$16000 (Tilt Elevate), \$19000 - \$20000 (Tilt Recline)

**Supplier:** Cubro

## Product Items on aSAP

Item Name	SIC	Supplier Code
Alltrack P3 12-16inWx12-18inD Tilt Elevate Powerchair	EM323	0600-0026
Alltrack P3 12-16inWx12-18inD Tilt Powerchair	EM322	0600-0024
Alltrack P3 12-16inWx12-18inD Tilt Recline Powerchair	EM324	0600-0025

## Product URL

<https://www.cubro.co.nz/healthcare-products/paediatric/wheelchairs/alltrack-p3-mid-wheel-drive-power-wheelchair>



# TA iQ Junior FWD Powerchair – Moderate Configuration



**Description:** TA iQ Junior Front Wheel Drive (FWD) power wheelchair. Base width: 600mm. Base length (without footplates): 820mm. Seat to floor height: 380mm. Turning radius: 450mm. Speed: 12.5km/h. Battery capacity: 80amps. Selections of: Tilt 45°, Tilt and elevate 11.8in, Tilt and recline with mechanical shear reduction.

**Cost bracket:** \$14000.00 - \$15000 (Tilt), \$16000.00 - \$17000 (Tilt Elevate), \$14000 - \$15000 (Tilt Recline)

**Supplier:** Cubro

**Intended Use:** Power wheelchair with moderate configuration (Paediatric). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1).

**Key Features:** Package includes: Standard controller and electronics (R-net expandable). Swing away controller mount / receiver. Centre mount leg hanger with flip up, angle adjustable footplates (2 individual pieces). Flip back, height adjustable arm supports with standard type, full length arm pads. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery charger. (Refer to each product package for specific components included)

**User Weight Capacity:** 140kg

# TA iQ Junior FWD Powerchair – Moderate Configuration



**Description:** TA iQ Junior Front Wheel Drive (FWD) power wheelchair. Base width: 600mm. Base length (without footplates): 820mm. Seat to floor height: 380mm. Turning radius: 450mm. Speed: 12.5km/h. Battery capacity: 80amps. Selections of: Tilt 45°, Tilt and elevate 11.8in, Tilt and recline with mechanical shear reduction.

**Cost bracket:** \$14000.00 - \$15000 (Tilt), \$16000.00 - \$17000 (Tilt Elevate), \$14000 - \$15000 (Tilt Recline)

**Supplier:** Cubro

## Product Items on aSAP

Item Name	SIC	Supplier Code
TA iQ FWD Jnr 10-22inWx10-24inD Tilt Elevate Powerchair	EM328	XTAIQFWDJ-TE
TA iQ FWD Jnr 10-22inWx10-24inD Tilt Powerchair	EM327	XTAIQFWDJ-T
TA iQ FWD Jnr 10-22inWx10-24inD Tilt Recline Powerchair	EM329	XTAIQFWDJ-TR

## Product URL

<https://www.cubro.co.nz/healthcare-products/powered-wheelchairs/front-wheel-drive/ta-iq-junior-wheelchair>

# Paediatric Power Wheelchair with High Configuration



Power Wheelchair Categories	Power Wheelchair with High Configuration
Paediatric	<p>Quantum Edge 3 Stretto (Allied Medical)</p> <p>Alltrack P3 (Cubro)</p> <p>TA IQ Junior (Cubro)</p>

**Green** = power wheelchairs in one category only

**Blue** = power wheelchairs in more than one category

MWD = Mid Wheel Drive

RWD = Rear Wheel Drive

FWD = Front Wheel Drive

# Stretto MWD Powerchair – High Configuration



**Description:** The narrowest, most compact and manoeuvrable power base in the Quantum range. 14in Drive Wheels. 5in Solid Castor Wheels Front. 6in Solid Castor Wheels Rear. Min Seat to Floor Height 17.5in. Maximum Rate Slope 6°. Max Speed 10km/h. Mid-Wheel Drive (MWD). Q-Logic 3 Expandable Electronics. SRS Suspension. High Efficiency And Mileage. 2 Pole Motor. Fender lights. Selections of: Power tilt, recline, elevate (TRE), power elevating centre mount (PCM) leg rest.

**Cost bracket:** > \$20000

**Supplier:** Allied Medical Ltd

**Intended Use:** Power wheelchair with high configuration (Paediatric, Adult). Suitable for indoor mobility and outdoor urban terrain. Designed for individuals with high and complex needs, including progressive conditions with the potential for changing needs. Features multiple power seat functions, supports alternative controls, and can support individuals displaying WMPM Level 2 clinical indicators.

**Key Features:** Package includes: Standard controller and electronic (Q-Logic 3 expandable). Swing away controller mount / receiver. Centre mount leg hanger with flip up, angle adjustable footplates. Flip back, height adjustable arm supports. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery and charger. Selections of: Power tilt, recline, elevate (TRE), power elevating centre mount (PCM) leg rest. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg

# Stretto MWD Powerchair – High Configuration



**Description:** The narrowest, most compact and manoeuvrable power base in the Quantum range. 14in Drive Wheels. 5in Solid Castor Wheels Front. 6in Solid Castor Wheels Rear. Min Seat to Floor Height 17.5in. Maximum Rate Slope 6°. Max Speed 10km/h. Mid-Wheel Drive (MWD). Q-Logic 3 Expandable Electronics. SRS Suspension. High Efficiency And Mileage. 2 Pole Motor. Fender lights. Selections of: Power tilt, recline, elevate (TRE), power elevating centre mount (PCM) leg rest.

**Cost bracket:** > \$20000

**Supplier:** Allied Medical Ltd

## Product Items on aSAP

Item Name	SIC	Supplier Code
Stretto MWD 16-20inWx12-16inD Tilt Recline (PCM) Powerchair	EM374	PKSTTRFPC16201216
Stretto MWD 16-20inWx16-20inD Tilt Recline (PCM) Powerchair	EM374	PKSTTRFPC16201620
Stretto MWD 16-20inWx20-22inD Tilt Recline (PCM) Powerchair	EM374	PKSTTRFPC16202024

Item Name	SIC	Supplier Code
Stretto MWD 12-16inWx12-16inD TRE Powerchair	EM373	PKSTTERFC12161216
Stretto MWD 14-16inWx16-20inD TRE (PCM) Powerchair	EM375	PKSTTERFPC14161620
Stretto MWD 14-16inWx16-20inD TRE Powerchair	EM373	PKSTTERFC14161620
Stretto MWD 16-20inWx16-20inD TRE (PCM) Powerchair	EM375	PKSTTERFPC16201620
Stretto MWD 16-20inWx16-20inD TRE Powerchair	EM373	PKSTTERFC16201620

## Product URL

<https://alliedmedical.co.nz/products/quantum-edge-3-stretto>

# Alltrack P3 MWD Powerchair – High Configuration



**Description:** Alltrack P3 Mid Wheel Drive (MWD) power wheelchair. Base width: 535mm. Base length: 890mm. Seat to floor height: 445mm. Turning radius: 510mm. Speed: 10km/h. Battery capacity: 50amps. Selections of: Power Tilt 50°, Power elevate 11in, Power recline with mechanical shear reduction seating functions, and Powered centre mount (PCM) elevating and articulating leg rest with flip-up, angle adjustable footplates.

**Cost bracket:** > \$20000

**Supplier:** Cubro

**Intended Use:** Power wheelchair with high configuration (Paediatric). Suitable for indoor mobility and outdoor urban terrain. Designed for individuals with complex needs and progressive conditions, with the potential for changing needs. Features multiple power seat functions, supports alternative controls, and can support individuals displaying WMPM Level 2 clinical indicators.

**Key Features:** Package includes: Standard controller and electronics (R-net expandable). Swing away controller mount / receiver. Centre mount leg hanger with flip up, angle adjustable footplates (2 individual pieces). Flip back; height adjustable arm supports with standard type, full-length arm pads. 2-point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery charger. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg

# Alltrack P3 MWD Powerchair – High Configuration



**Description:** Alltrack P3 Mid Wheel Drive (MWD) power wheelchair. Base width: 535mm. Base length: 890mm. Seat to floor height: 445mm. Turning radius: 510mm. Speed: 10km/h. Battery capacity: 50amps. Selections of: Power Tilt 50°, Power elevate 11in, Power recline with mechanical shear reduction seating functions, and Powered centre mount (PCM) elevating and articulating leg rest with flip-up, angle adjustable footplates.

**Cost bracket:** > \$20000

**Supplier:** Cubro

## Product Items on aSAP

Item Name	SIC	Supplier Code
Alltrack P3 12-16inWx12-18inD TRE (PCM) Powerchair	EM326	0600-0028
Alltrack P3 12-16inWx12-18inD TRE Powerchair	EM325	0600-0027

## Product URL

<https://www.cubro.co.nz/healthcare-products/paediatric/wheelchairs/alltrack-p3-mid-wheel-drive-power-wheelchair>



# TA iQ Junior FWD Powerchair – High Configuration



**Description:** TA iQ Junior Front Wheel Drive (FWD) power wheelchair. Base width: 600mm. Base length (without footplates): 820mm. Seat to floor height: 380mm. Turning radius: 450mm. Speed: 12.5km/h. Battery capacity: 80amps. Selections of: Power tilt 45°, Power elevate 11.8in, Power recline with mechanical shear reduction, and Powered centre mount (PCM) elevating and articulating leg rest with flip-up, angle adjustable footplates.

**Cost bracket:** \$19000 - \$20000

**Supplier:** Cubro

**Intended Use:** Power wheelchair with high configuration (Paediatric). Suitable for indoor mobility and outdoor urban terrain. Designed for individuals with high and complex needs, including progressive conditions with the potential for changing needs. Features multiple power seat functions, supports alternative controls, and can support individuals displaying WMPM Level 2 clinical indicators.

**Key Features:** Package includes: Standard controller and electronics (R-net expandable). Swing away controller mount / receiver. Centre mount leg hanger with flip up, angle adjustable footplates (2 individual pieces). Flip back, height adjustable arm supports with standard type, full length arm pads. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery charger. (Refer to each product package for specific components included)

**User Weight Capacity:** 140kg



# TA iQ Junior FWD Powerchair – High Configuration



**Description:** TA iQ Junior Front Wheel Drive (FWD) power wheelchair. Base width: 600mm. Base length (without footplates): 820mm. Seat to floor height: 380mm. Turning radius: 450mm. Speed: 12.5km/h. Battery capacity: 80amps. Selections of: Power tilt 45°, Power elevate 11.8in, Power recline with mechanical shear reduction, and Powered centre mount (PCM) elevating and articulating leg rest with flip-up, angle adjustable footplates.

**Cost bracket:** \$19000 - \$20000

**Supplier:** Cubro

## Product Items on aSAP

Item Name	SIC	Supplier Code
TA iQ FWD Jnr 10-22inWx10-24inD TRE (PCM) Powerchair	EM331	XTAIQFWDJ-MF
TA iQ FWD Jnr 10-22inWx10-24inD TRE Powerchair	EM330	XTAIQFWDJ-TRE

## Product URL

<https://www.cubro.co.nz/healthcare-products/powered-wheelchairs/front-wheel-drive/ta-iq-junior-wheelchair>

# Accessible List

## Adult Power Wheelchair Range



# Adult Power Wheelchair Products List Overview



Power Wheelchair Categories	Power Wheelchair with Minimal Configuration	Power Wheelchair with Moderate Configuration	Power Wheelchair with High Configuration
Adult	<p><u>MWD options:</u> M1 (non-tilt) (Permobil)</p> <p><u>RWD options:</u> Karma KP25 2.0 (Allied Medical) Q100 RWD (Medifab) Q200 RWD (Medifab)</p>	<p><u>MWD options:</u> M1 (with tilt) (Permobil) Pride J623 2.0 (Allied Medical) Q500M (Medifab) TA IQ MWD (Cubro) TDXSP2 REM216 (Invacare) Quantum Edge 3 (Allied Medical) M3 Corpus (Permobil)</p> <p>Q500 Mini (Medifab) Quantum Edge 3 Stretto (Allied Medical)</p> <p><u>RWD options:</u> Aviva RX REM216 (Invacare)</p>	<p><u>MWD options:</u> Q500M (Medifab) TA IQ MWD (Cubro) TDXSP2 REM400 (Invacare) Quantum Edge 3 (Allied Medical) M3 Corpus (Permobil) Q500 Mini (Medifab) Quantum Edge 3 Stretto (Allied Medical)</p> <p><u>RWD options:</u> Aviva RX REM400 (Invacare)</p> <p><u>FWD options:</u> F3 Corpus (Permobil)</p>

**Green** = power wheelchairs in one category only

**Blue** = power wheelchairs in more than one category

MWD = Mid Wheel Drive

RWD = Rear Wheel Drive

FWD = Front Wheel Drive

# Adult Power Wheelchair Range

Accessable has categorised adult power wheelchair range according to the level of configuration and drive wheel types.

Level of Configuration
Minimal configuration
Moderate configuration
High configuration

Drive Wheel Type
Mid Wheel Drive (MWD)
Rear Wheel Drive (RWD)
Front Wheel Drive (MWD)

# Adult Power Wheelchair Range

Within this range, there are power wheelchair options with:

- **Minimal configuration through to highly configurable options** that meet clients' needs, ranging from general use to high and complex requirements
- **A compact base and a small turning radius** for use in indoor environments with confined spaces
- **Seat functions**, ranging from non-tilt options to multiple power seat functions
- **Controller capabilities**, non-expandable/ expandable electronics
- **Varying degrees of motor/ battery capacities**

*Note: When selecting power wheelchairs and features, it's important to consider funding criteria. Assessors should determine the essential needs based on the client's context, such as the necessity for a power wheelchair with extended motor capacity for community mobility.*

# Adult Power Wheelchair Range

## Differences between power wheelchair front, mid and rear drive wheel type

- Useful summary document: [Comparison of Front, Mid, and Rear Wheel Drive Power Chairs](#)



- **Website link:**
  - [https://www.health.qld.gov.au/\\_data/assets/pdf\\_file/0028/428482/pdw-c-comparison.pdf](https://www.health.qld.gov.au/_data/assets/pdf_file/0028/428482/pdw-c-comparison.pdf)

# Adult Power Wheelchair with Minimal Configuration



Power Wheelchair Categories	Power Wheelchair with Minimal Configuration
Adult	<p><u>MWD options:</u> M1 (non-tilt) (Permobil)</p> <p><u>RWD options:</u> Karma KP25 2.0 (Allied Medical) Q100 RWD (Medifab) Q200 RWD (Medifab)</p>

Green = power wheelchairs in one category only

Blue = power wheelchairs in more than one category

MWD = Mid Wheel Drive

RWD = Rear Wheel Drive

FWD = Front Wheel Drive

# Karma KP25.2 Low Powerchair



**Description:** KP25.2 Powerchair 19In  
Low STF

**Cost bracket:** **\$5000.00 - \$6000**

**Supplier:** Allied Medical Ltd

**Intended Use:** Power Wheelchair with low configuration (Adult). Suitable for indoor mobility and immediate outdoor environment around home. This is for individuals new to using power wheelchairs and can support individuals displaying WMPM Level 1 indicators.

**Key Features:** Package includes: Battery charger, batteries and positioning belt. (Refer to each product package for specific components included)

**User Weight Capacity:** 135kg



# Karma KP25.2 Low Powerchair



**Description:** KP25.2 Powerchair 19In Low STF

**Cost bracket:** \$5000.00 - \$6000

**Supplier:** Allied Medical Ltd

## Product Items on aSAP

Item Name	SIC	Supplier Code
Karma KP25.2 Low 14inx16in Powerchair	EM103	AKP2521416L
Karma KP25.2 Low 16inx16in Powerchair	EM103	AKP2521616L
Karma KP25.2 Low 16inx17in Powerchair	EM103	KP2521616L
Karma KP25.2 Low 16inx18in Powerchair	EM103	AKP2521618L
Karma KP25.2 Low 18inx16in Powerchair	EM103	AKP2521816L
Karma KP25.2 Low 18inx18in Powerchair	EM103	AKP2521818L
Karma KP25.2 Low 20inx18in Powerchair	EM103	AKP2522018L

## Product URL

<https://alliedmedical.co.nz/products/karma-kp252-lightweight-folding-powerchair>

# Permobil M1 (Non-Tilt) Powerchair



**Description:** Permobil M1 Mid-Wheel-Drive Power Wheelchair (Non-Tilt). Max Speed 8km/h, 14in Drive Wheel, Floor-to-Seat Height - 16, 17 or 18in.

**Cost bracket:** \$6000.00 - \$7000

**Supplier:** Permobil New Zealand Limited

**Intended Use:** Power Wheelchair with low configuration (Adult). Suitable for indoor mobility and immediate outdoor environment around home. This is for individuals new to using power wheelchairs and can support individuals displaying WMPM Level 1 indicators.

**Key Features:** Package includes: M1 Power Base, Seat Pan only (no power tilt). Must choose seat width size width. Depth is fully adjustable 14-22in. R-net VR2 Non-expandable Joystick, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Locking Flip Back Height Adjustable Armrests, Manual Elevating Centre-Mount Legrest with Calf Supports. Tubular Backrest or Curved Solid Backrest, Standard 2-point Positioning Belt, Transfer Handles, LED Light Kit, Leatherette or Fabric Armrests. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg

# Permobil M1 (Non-Tilt) Powerchair



**Description:** Permobil M1 Mid-Wheel-Drive Power Wheelchair (Non-Tilt). Max Speed 8km/h, 14in Drive Wheel, Floor-to-Seat Height - 16, 17 or 18in.

**Cost bracket:** **\$6000.00 - \$7000**

**Supplier:** Permobil New Zealand Limited

## Product Items on aSAP

Item Name	SIC	Supplier Code
Permobil M1 14inWx14-22inD (Non-Tilt) Powerchair	EM300	PKGM1
Permobil M1 16inWx14-22inD (Non-Tilt) Powerchair	EM301	PKGM1
Permobil M1 18inWx14-22inD (Non-Tilt) Powerchair	EM302	PKGM1
Permobil M1 20inWx14-22inD (Non-Tilt) Powerchair	EM303	PKGM1
Permobil M1 22inWx14-22inD (Non-Tilt) Powerchair	EM304	PKGM1

**Note:** Seat width is fixed and cannot be adjusted. If a client's seat width changes, a new seat top of the required width can be installed.

## Product URL

<https://www.permobil.com/en-nz/products/power-wheelchairs/m1>

# Q100R Powerchair



**Description:** Quickie Q100R Standard RWD Powerchair, No Power Seat Functions, Includes Seat, Back Upholstery, Armrests & Footrests. 12in Drive Wheels. 8in Castors. 6KM Speed. 31KM Range.

**Cost bracket:** \$4000.00 - \$5000

**Supplier:** Medifab Limited

**Intended Use:** Power Wheelchair with low configuration (Adult). Suitable for indoor mobility and immediate outdoor environment around home. This is for individuals new to using power wheelchairs and can support individuals displaying WMPM Level 1 indicators.

**Key Features:** Package includes: VR2 Non-Expandable Electronics. Standard Controller. Swing Away Controller Mount. 12in Drive Wheels. 8in Castors. Transit Tie Downs. 2 Pole Motors - 6KM. 55Ah Gel Batteries. Battery Charger. Removable Armrests. Swing Away Legrests with Composite Footplates. Heel Loops. Padded Seat Cushion. Back Upholstery. Standard Lap Belt. Seat width 17-19.5in (adjusted through arm rests)

**User Weight Capacity:** 125kg

## Product Items on aSAP

Item Name	SIC	Supplier Code
Q100R 18inWx16-20inD Powerchair	EM337	1636-1001-000

## Product URL

<https://medifab.com/nz/product/quickie-q100-r-powered-wheelchair/>

# Q200R Powerchair



**Description:** Quickie Q200R Standard RWD Powerchair, No Power Seat Functions, Includes Seat, Back Upholstery, Armrests & Footrests. 14in Drive Wheels. 9in Castors. 10KM Speed. 31KM Range.

**Cost bracket:** **\$5000.00 - \$6000**

**Supplier:** Medifab Limited

**Intended Use:** Power Wheelchair with low configuration (Adult). Suitable for indoor mobility and immediate outdoor environment around home. This is for individuals new to using power wheelchairs and can support individuals displaying WMPM Level 1 indicators.

**Key Features:** Package includes: VR2 Non-Expandable Electronics. Standard Controller. Swing Away Controller Mount. 14in Drive Wheels. 9in Castors. Transit Tie Downs. 2 Pole Motors - 10KM. 55Ah Gel Batteries. Battery Charger. Removable Armrests. Swing Away Legrests with Composite Footplates. Heel Loops. Padded Seat Cushion. Back Upholstery. Standard Lap Belt. Seat width 17-19.5in (adjusted through arm rests)

**User Weight Capacity:** 136kg

## Product Items on aSAP

Item Name	SIC	Supplier Code
Q200R 18inWx16-20inD Powerchair	EM338	1637-1001-000

## Product URL

<https://medifab.com/product/quickie-q200-r-powered-wheelchair/>

# Adult Power Wheelchair with Moderate Configuration



Power Wheelchair Categories	Power Wheelchair with Moderate Configuration
Adult	<p><b><u>MWD options:</u></b></p> <p>M1 (with tilt) (Permobil)</p> <p>Pride J623 2.0 (Allied Medical)</p> <p>Q500M (Medifab)</p> <p>TA IQ MWD (Cubro)</p> <p>TDXSP2 REM216 (Invacare)</p> <p>Quantum Edge 3 (Allied Medical)</p> <p>M3 Corpus (Permobil)</p> <p>Q500 Mini (Medifab)</p> <p>Quantum Edge 3 Stretto (Allied Medical)</p> <p><b><u>RWD options:</u></b></p> <p>Aviva RX REM216 (Invacare)</p>

Green = power wheelchairs in one category only

Blue = power wheelchairs in more than one category

MWD = Mid Wheel Drive

RWD = Rear Wheel Drive

FWD = Front Wheel Drive



# Permobil M1 Tilt Powerchair



**Description:** Permobil M1 Mid-Wheel-Drive Power Wheelchair Tilt. Max Speed 8km/h, 14in Drive Wheel, Floor-to-Seat Height - 16, 17 or 18in. 50° Power Tilt.

**Cost bracket:** **\$11000.00 - \$12000**

**Supplier:** Permobil New Zealand Limited

**Intended Use:** Power wheelchair with moderate configuration (Adult). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1).

**Key Features:** Package includes: M1 Power Base with 50° Power Tilt, Seat Pan, must choose seat width size width. Depth is fully adjustable 14-22in. R-net LED Remote Non-expandable Joystick, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Locking Flip Back Height Adjustable Armrests, Manual Elevating Centre-Mount Legrest with Calf Supports. Tubular Backrest or Curved Solid Backrest, Standard 2-point Positioning Belt, Transfer Handles, LED Light Kit, Leatherette or Fabric Armrests. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg



# Permobil M1 Tilt Powerchair



**Description:** Permobil M1 Mid-Wheel-Drive Power Wheelchair Tilt. Max Speed 8km/h, 14in Drive Wheel, Floor-to-Seat Height - 16, 17 or 18in. 50° Power Tilt.

**Cost bracket:** **\$11000.00 - \$12000**

**Supplier:** Permobil New Zealand Limited

## Product Items on aSAP

Item Name	SIC	Supplier Code
Permobil M1 14inWx14-22inD Tilt Powerchair	EM305	PKGM1T
Permobil M1 16inWx14-22inD Tilt Powerchair	EM306	PKGM1T
Permobil M1 18inWx14-22inD Tilt Powerchair	EM307	PKGM1T
Permobil M1 20inWx14-22inD Tilt Powerchair	EM308	PKGM1T
Permobil M1 22inWx14-22inD Tilt Powerchair	EM309	PKGM1T

**Note:** Seat width is fixed and cannot be adjusted. If a client's seat width changes, a new seat top of the required width can be installed.

## Product URL

<https://www.permobil.com/en-nz/products/power-wheelchairs/m1>

# Jazzy J623 2.0 MWD Powerchair



**Description:** Ideal for users who like to be active with enhanced suspension for a comfortable ride. Easy to use and cost-effective. Mid-Wheel Drive (MWD). VR2 Electronics. Active-Trak ATX Adjustable Suspension. High Efficiency And Mileage. Linux 2 Pole Motor. Quick Height. Adjustable Armrests.

**Cost bracket:** \$15000.00 - \$16000 (Tilt), \$17000.00 - \$18000 (Tilt Elevate)

**Supplier:** Allied Medical Ltd

**Intended Use:** Power wheelchair with moderate configuration (Adult). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1).

**Key Features:** Package consists of: Standard controller and electronic (VR2 non expandable). Swing away controller mount / receiver. Centre mount leg hanger with flip up, angle adjustable Footplates. Flip back, height adjustable arm supports. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery and charger. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg

# Jazzy J623 2.0 MWD Powerchair



**Description:** Ideal for users who like to be active with enhanced suspension for a comfortable ride. Easy to use and cost-effective. Mid-Wheel Drive (MWD). VR2 Electronics. Active-Trak ATX Adjustable Suspension. High Efficiency And Mileage. Linux 2 Pole Motor. Quick Height. Adjustable Armrests.

**Cost bracket:** \$15000.00 - \$16000 (Tilt), \$17000.00 - \$18000 (Tilt Elevate)

**Supplier:** Allied Medical Ltd

## Product Items on aSAP

Item Name	SIC	Supplier Code
Jazzy J623 2.0 MWD 12-16inWx12-16inD Tilt Powerchair	EM388	PKJ6232TFC12161216
Jazzy J623 2.0 MWD 12-16inWx16-20inD Tilt Powerchair	EM388	PKJ6232TFC12161620
Jazzy J623 2.0 MWD 16-20inWx16-20inD Tilt Powerchair	EM388	PKJ6232TFC16201620
Jazzy J623 2.0 MWD 16-20inWx20-22inD Tilt Powerchair	EM388	PKJ6232TFC16202024
Jazzy J623 2.0 MWD 20-24inWx16-20inD Tilt Powerchair	EM388	PKJ6232TFC20241620
Jazzy J623 2.0 MWD 20-24inWx20-22inD Tilt Powerchair	EM388	PKJ6232TFC20242024

Item Name	SIC	Supplier Code
Jazzy J623 2.0 MWD 14-16inWx16-20inD Tilt Elevate Powerchair	EM419	PKJ6232TEFC14161620
Jazzy J623 2.0 MWD 16-20inWx16-20inD Tilt Elevate Powerchair	EM420	PKJ6232TEFC16201620
Jazzy J623 2.0 MWD 16-20inWx20-24inD Tilt Elevate Powerchair	EM421	PKJ6232TEFC16202024
Jazzy J623 2.0 MWD 20-24inWx16-20inD Tilt Elevate Powerchair	EM422	PKJ6232TEFC20241620
Jazzy J623 2.0 MWD 20-24inWx20-24inD Tilt Elevate Powerchair	EM423	PKJ6232TEFC20242024

## Product URL

<https://alliedmedical.co.nz/products/jazzy-623-20>

# Q500M MWD Powerchair – Moderate Configuration



**Description:** Q500M MWD Powerchair. Includes Adjustable Seat, Armrests & Footrests. 13in Drive Wheels. 7in Castors. 60 Ah Batteries. 10km Speed. 40km Range. Options: 50 Degree Posterior Tilt. 30cm Seat Elevate. 130 Degree Recline. Swing Away Angle Adjustable Footplates or Centre-mount (CM) Footrests.

**Cost bracket:** \$11000.00 - \$12000 (Tilt), \$14000.00 - \$15000 (Tilt Elevate), \$15000 - \$16000 (Tilt Recline)

**Supplier:** Medifab Limited

**Intended Use:** Power wheelchair with moderate configuration (Adult). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1).

**Key Features:** Package includes: R-Net Expandable Electronics. 13in Drive Wheels. 7in Castors. Transit Tie Downs. 4 Pole Motors. 60Ah Gel Batteries. Battery Charger. Sedeo Lite Seat Top. Power Seat Tilt. 70 Degree Swing Away Legrests. Aluminium Angle Adjustable Footplates. Or Centre-mount (CM) Footrests. Cantilever Flip Back Armrests. RHS Swing Away Controller Mount. CJSM1 Standard Joystick. Standard Lap Belt. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg

# Q500M MWD Powerchair – Moderate Configuration



**Description:** Q500M MWD Powerchair. Includes Adjustable Seat, Armrests & Footrests. 13in Drive Wheels. 7in Castors. 60 Ah Batteries. 10km Speed. 40km Range. Options: 50 Degree Posterior Tilt. 30cm Seat Elevate. 130 Degree Recline. Swing Away Angle Adjustable Footplates or Centre-mount (CM) Footrests.

**Cost bracket:** \$11000.00 - \$12000 (Tilt), \$14000.00 - \$15000 (Tilt Elevate), \$15000 - \$16000 (Tilt Recline)

**Supplier:** Medifab Limited

## Product Items on aSAP

Item Name	SIC	Supplier Code
Q500M MWD 16-20inWx16-22inD Tilt (CM) Powerchair	EM343	1652-1901-004
Q500M MWD 16-20inWx16-22inD Tilt (SA) Powerchair	EM339	1652-1901-001
Item Name	SIC	Supplier Code
Q500M MWD 16-20inWx16-22inD Tilt Recline (CM) Powerchair	EM344	1652-1901-005
Q500M MWD 16-20inWx16-22inD Tilt Recline (SA) Powerchair	EM340	1652-1901-002
Item Name	SIC	Supplier Code
Q500M MWD 16-20inW-16-22inD Tilt Elevate (CM) Powerchair	EM345	1652-1901-006
Q500M MWD 16-20inW-16-22inD Tilt Elevate (SA) Powerchair	EM341	1652-1901-003

## Product URL

<https://medifab.com/nz/product/quickie-q500m/>

# TA iQ MWD Powerchair – Moderate Configuration



**Description:** TA iQ Mid Wheel Drive (MWD) power wheelchair. Base width: 630mm. Base length (without footplates): 910mm. Seat to floor height: 380mm. Turning radius: 450mm. Speed: 12.5km/h. Battery capacity: 80amps. Selections of: Tilt, Tilt and elevate, Tilt and recline seating functions

**Cost bracket:** \$14000.00 - \$15000 (Tilt), \$16000.00 - \$17000 (Tilt Elevate), \$16000 - \$17000 (Tilt Recline)  
**Supplier:** Cubro

**Intended Use:** Power wheelchair with moderate configuration (Adult). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1).

**Key Features:** Package includes: Standard controller and electronics (R-net expandable). Swing away controller mount / receiver. Centre mount leg hanger with flip up, angle adjustable footplates (2 individual pieces). Flip back, height adjustable arm supports with standard type, full length arm pads. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery charger. (Refer to each product package for specific components included)

**User Weight Capacity:** 140kg



# TA iQ MWD Powerchair – Moderate Configuration



**Description:** TA iQ Mid Wheel Drive (MWD) power wheelchair. Base width: 630mm. Base length (without footplates): 910mm. Seat to floor height: 380mm. Turning radius: 450mm. Speed: 12.5km/h. Battery capacity: 80amps. Selections of: Tilt, Tilt and elevate, Tilt and recline seating functions

**Cost bracket:** \$14000.00 - \$15000 (Tilt), \$16000.00 - \$17000 (Tilt Elevate), \$16000 - \$17000 (Tilt Recline)  
**Supplier:** Cubro

## Product Items on aSAP

Item Name	SIC	Supplier Code
TA iQ MWD 16-24inWx10-24inD Tilt Elevate Powerchair	EM333	XTAIQMWD-TE
TA iQ MWD 16-24inWx10-24inD Tilt Powerchair	EM332	XTAIQMWD-T
TA iQ MWD 16-24inWx10-24inD Tilt Recline Powerchair	EM334	XTAIQMWD-TR

## Product URL

<https://www.cubro.co.nz/healthcare-products/powered-wheelchairs/mid-wheel-drive/ta-iq-mwd-wheelchair>

# TDXSP2 MWD Powerchair – Moderate Configuration



**Description:** TDXSP2 Narrow or Wide Mid wheel drive (MWD) base, Selections of: Tilt, Tilt and elevate, Tilt and recline seating functions. Non-expandable joystick. Optional Centre-mount or 70 Degree Swing away leg rests.

**Cost bracket:** \$13000.00 - \$14000 (Tilt), \$18000.00 - \$19000 (Tilt Elevate or Tilt Recline)

**Supplier:** Invacare New Zealand

**Intended Use:** Power wheelchair with moderate configuration (Adult). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1).

**Key Features:** Package includes: TDXSP2 Narrow or Wide Base, ULM seating system (selections of: Tilt, Tilt and elevate, Tilt and recline), 2-point padded hip belt, REM216 non-expandable Joystick, 52amp batteries, Centre-mount Leg rest or 70 Degree swing away leg rests, fixed mount joystick, Dual-Post Flip-Back Armrests (Non-Recline Seats) or Cantilever Armrests (Recline seat). (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg



# TDXSP2 MWD Powerchair – Moderate Configuration



**Description:** TDXSP2 Narrow or Wide Mid wheel drive (MWD) base, Selections of: Tilt, Tilt and elevate, Tilt and recline seating functions. Non-expandable joystick. Optional Centre-mount or 70 Degree Swing away leg rests.

**Cost bracket:** \$13000.00 - \$14000 (Tilt), \$18000.00 - \$19000 (Tilt Elevate or Tilt Recline)

**Supplier:** Invacare New Zealand

## Product Items on aSAP

Item Name	SIC	Supplier Code
TDXSP2 MWD Nrrw 14-24inWx14-24inD Tilt Elevate Powerchair	EM409	PKG-TDX NARR 3
TDXSP2 MWD Nrrw 14-24inWx14-24inD Tilt Powerchair	EM407	PKG-TDX NARR 1
TDXSP2 MWD Nrrw 14-24inWx14-24inD Tilt Recline Powerchair	EM408	PKG-TDX NARR 2
Item Name	SIC	Supplier Code
TDXSP2 MWD Wide 14-24inWx14-24inD Tilt Elevate Powerchair	EM412	PKG-TDX WIDE 3
TDXSP2 MWD Wide 14-24inWx14-24inD Tilt Powerchair	EM410	PKG-TDX WIDE 1
TDXSP2 MWD Wide 14-24inWx14-24inD Tilt Recline Powerchair	EM411	PKG-TDX WIDE 2

## Product URL

<https://www.invacare.co.nz/move/powerd-mobility/powerd-wheelchairs/tdx-sp2-linx>

# Edge 3 MWD Powerchair – Moderate Configuration



**Description:** A Mid Wheel Drive (MWD) powerchair, easy to customise the mechanical and electronic set up. Ideal for home and in the community. 14in Drive Wheels. 6in Solid Castor Wheels Front. 6in Solid Castor Wheels Rear. Min Seat to Floor Height 16.5in. Maximum Rate Slope 7.5°. Max Speed 10km/h. Selections of: Power tilt. Power tilt and elevate. Power Tilt and Recline.

**Cost bracket:** \$16000.00 - \$17000 (Tilt), \$19000.00 - \$20000 (Tilt Elevate), >\$20000 (Tilt Recline)

**Supplier:** Allied Medical Ltd

**Intended Use:** Power wheelchair with moderate configuration (Adult). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1).

**Key Features:** Package includes: Standard controller and electronic (Q-Logic 3 expandable. Swing away controller mount / receiver. Centre mount leg hanger with flip up, angle adjustable footplates. Flip back, height adjustable arm supports. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery and charger. Selections of: Power tilt. Power tilt and elevate. Power Tilt and Recline. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg

# Edge 3 MWD Powerchair – Moderate Configuration



## Product Items on aSAP

Item Name	SIC	Supplier Code
Edge 3 MWD 12-16inWx16-20inD Tilt Powerchair	EM376	PKE3TFC12161620
Edge 3 MWD 16-20inWx12-16inD Tilt Powerchair	EM376	PKE3TFC16201216
Edge 3 MWD 16-20inWx20-22inD Tilt Powerchair	EM376	PKE3TFC16202024
Edge 3 MWD 20-22inWx16-20inD Tilt Powerchair	EM376	PKE3TFC20241620
Edge 3 MWD 20-22inWx20-22inD Tilt Powerchair	EM376	PKE3TFC20242024

## Product Items on aSAP

Item Name	SIC	Supplier Code
Edge 3 MWD 14-16inWx16-20inD Tilt Elevate Powerchair	EM377	PKE3TFC20242024
Edge 3 MWD 16-20inWx20-22inD Tilt Elevate Powerchair	EM377	PKE3TEFC16202024
Edge 3 MWD 20-22inWx20-22inD Tilt Elevate Powerchair	EM377	AEDGE3REAC2024A
Edge 3 MWD 20-22inWx16-20inD Tilt Elevate Powerchair	EM377	PKE3TEFC20241620

Item Name	SIC	Supplier Code
Edge 3 MWD 12-16inWx12-16inD Tilt Recline Powerchair	EM378	PKE3TRFC12161216
Edge 3 MWD 12-16inWx16-20inD Tilt Recline Powerchair	EM378	PKE3TRFC12161620
Edge 3 MWD 16-20inWx12-16inD Tilt Recline Powerchair	EM378	PKE3TRFC16201216
Edge 3 MWD 16-20inWx20-22inD Tilt Recline Powerchair	EM378	PKE3TRFC16202024
Edge 3 MWD 20-22inWx16-20inD Tilt Recline Powerchair	EM378	PKE3TRFC20241620
Edge 3 MWD 20-22inWx20-22inD Tilt Recline Powerchair	EM378	PKE3TRFC20242024

## Product URL

<https://alliedmedical.co.nz/products/quantum-edge-3>

# Permobil M3 Powerchair - Moderate Configuration



**Description:** Permobil M3 Mid-Wheel-Drive Power Wheelchair with Power Tilt and Power Elevate. Max Speed 10km/h, 14in Drive Wheel, Floor-to-Seat Height - 17.5-19.5in. Selections of: 50° Power Tilt only, Power Tilt and 12in Power Elevate, or Power Tilt and 180° Power Recline.

**Cost bracket:** \$16000.00 - \$17000 (Tilt), \$19000.00 - \$20000 (Tilt Elevate or Tilt Recline)

**Supplier:** Permobil New Zealand Limited

**Intended Use:** Power wheelchair with moderate configuration (Adult). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1).

**Key Features:** Package includes: M3 Power Base. Selections of: 50° Power Tilt only, Power Tilt and 12in Power Elevate, or Power Tilt and 180° Power Recline, Seat Pan, Width adjustable 17-23in, Depth adjustable 14-22in. R-net Colour CJSM Expandable Joystick, ICS Seating Switchbox, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Anti-Shear Backrest Assembly, Flip Back Height and Width Adjustable Armrests, Manual Elevating Centre-Mount Leg rest, with Calf Supports, USB Device Charger. Standard 2-point Positioning Belt, Transfer Handles, LED Light Kit, Leatherette or Fabric Armrests. (Refer to each product package for specific components included)

**User Weight Capacity:** 150kg

# Permobil M3 Powerchair - Moderate Configuration



**Description:** Permobil M3 Mid-Wheel-Drive Power Wheelchair with Power Tilt and Power Elevate. Max Speed 10km/h, 14in Drive Wheel, Floor-to-Seat Height - 17.5-19.5in. Selections of: 50° Power Tilt only, Power Tilt and 12in Power Elevate, or Power Tilt and 180° Power Recline.

**Cost bracket:** \$16000.00 - \$17000 (Tilt), \$19000.00 - \$20000 (Tilt Elevate or Tilt Recline)

**Supplier:** Permobil New Zealand Limited

## Product Items on aSAP

Item Name	SIC	Supplier Code
Permobil M3 17-23inWx14-22inD Tilt Powerchair	EM310	PKGM3T
Item Name	SIC	Supplier Code
Permobil M3 17-23inWx14-22inD Tilt Elevate Powerchair	EM312	PKGM3TE
Permobil M3 17-23inWx14-22inD Tilt Recline Powerchair	EM314	PKGM3TR

## Product URL

<https://www.permobil.com/en-nz/products/power-wheelchairs/m3-corpus>

# Q500 Mini MWD Powerchair – Moderate Configuration



**Description:** Q500 Mini MWD Powerchair. 13in Drive Wheels. 7in Castors. 50 Ah Batteries. 10km Speed. 30km Range. Selections of: 50 Degree Posterior Tilt. Tilt and 30cm Seat Elevate. Tilt and 130 Degree Recline. Single Flip Up Centre-mount (CM) Footrests or Manual Swing Away (SA) Footrests with Angle Adjustable Footplates.

**Cost bracket:** \$13000.00 - \$15000 (Tilt), \$15000.00 - \$16000 (Tilt Elevate), \$17000 - \$19000 (Tilt Recline)

**Supplier:** Medifab Limited

**Intended Use:** Power wheelchair with moderate configuration (Adult). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1).

**Key Features:** Package includes: R-Net Expandable Electronics. 13in Drive Wheels. 7in Castors. Transit Tie Downs. 4 Pole Motors. 50Ah Gel Batteries. Battery Charger. Selections of: Tilt only. Tilt and Elevate. Tilt and Recline. Sedeo Lite Seat Top (Non-recline) or Sedeo Pro Seat Top (Power Recline). Single Flip Up Centre-mount (CM) Footrests or 70 Degree Swing Away Legrest with Aluminium Angle Adjustable Footplates. Cantilever Flip Back Armrests. RHS Swing Away Controller Mount. CJSM1 Standard Joystick. Standard Lap Belt. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg



# Q500 Mini MWD Powerchair – Moderate Configuration



**Description:** Q500 Mini MWD Powerchair. 13in Drive Wheels. 7in Castors. 50 Ah Batteries. 10km Speed. 30km Range. Selections of: 50 Degree Posterior Tilt. Tilt and 30cm Seat Elevate. Tilt and 130 Degree Recline. Single Flip Up Centre-mount (CM) Footrests or Manual Swing Away (SA) Footrests with Angle Adjustable Footplates.

**Cost bracket:** **\$13000.00 - \$15000 (Tilt), \$15000.00 - \$16000 (Tilt Elevate), \$17000 - \$19000 (Tilt Recline)**

**Supplier:** Medifab Limited

## Product Items on aSAP

Item Name	SIC	Supplier Code
Q500 Mini MWD 16-20inWx16-22inD Tilt (CM) Powerchair	EM354	1651-1901-004
Q500 Mini MWD 16-20inWx16-22inD Tilt (SA) Powerchair	EM351	1651-1901-001
Q500 Mini MWD 16-20inWx16-22inD Tilt Elevate (CM) Powerchair	EM356	1651-1901-006
Q500 Mini MWD 16-20inWx16-22inD Tilt Elevate (SA) Powerchair	EM353	1651-1901-003
Q500 Mini MWD 16-20inWx16-22inD Tilt Recline (CM) Powerchair	EM355	1651-1901-005
Q500 Mini MWD 16-20inWx16-22inD Tilt Recline (SA) Powerchair	EM352	1651-1901-002

## Product URL

<https://medifab.com/nz/product/quickie-q500m-mini/>

# Stretto MWD Powerchair – Moderate Configuration



**Description:** The narrowest, most compact and manoeuvrable power base in the Quantum range. 14in Drive Wheels. 5in Solid Castor Wheels Front. 6in Solid Castor Wheels Rear. Min Seat to Floor Height 17.5in. Maximum Rate Slope 6°. Max Speed 10km/h. Mid-Wheel Drive (MWD). Q-Logic 3 Expandable Electronics. SRS Suspension. High Efficiency And Mileage. 2 Pole Motor. Fender lights. Selections of: Power tilt. Power tilt and elevate. Power Tilt and Recline.

**Cost bracket:** \$15000.00 - \$16000 (Tilt), \$18000.00 - \$19000 (Tilt Elevate), \$19000.00 - \$20000 (Tilt Recline)

**Supplier:** Allied Medical Ltd

**Intended Use:** Power wheelchair with moderate configuration (Paediatric, Adult). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1).

**Key Features:** Package includes: Standard controller and electronic (Q-Logic 3 expandable). Swing away controller mount / receiver. Centre mount leg hanger with flip up, angle adjustable footplates. Flip back, height adjustable arm supports. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery and charger. Selections of: Power tilt. Power tilt and elevate. Power Tilt and Recline. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg



# Stretto MWD Powerchair – Moderate Configuration



## Product Items on aSAP

Item Name	SIC	Supplier Code
Stretto MWD 12-16inWx12-16inD Tilt Powerchair	EM370	PKSTTFC12161216
Stretto MWD 12-16inWx16-20inD Tilt Powerchair	EM370	PKSTTFC12161620
Stretto MWD 16-20inWx12-16inD Tilt Powerchair	EM370	PKSTTFC16201216
Stretto MWD 16-20inWx16-20inD Tilt Powerchair	EM370	PKSTTFC16201620
Stretto MWD 16-20inWx20-22inD Tilt Powerchair	EM370	PKSTTFC16202024

## Product Items on aSAP

Item Name	SIC	Supplier Code
Stretto MWD 12-16inWx12-16inD Tilt Elevate Powerchair	EM371	PKSTTEFC12161216
Stretto MWD 14-16inWx16-20inD Tilt Elevate Powerchair	EM371	PKSTTEFC14161620
Stretto MWD 16-20inWx16-20inD Tilt Elevate Powerchair	EM371	PKSTTEFC16201620

Item Name	SIC	Supplier Code
Stretto MWD 12-16inWx12-16inD Tilt Recline Powerchair	EM372	PKSTTRFC12161216
Stretto MWD 12-16inWx16-20inD Tilt Recline Powerchair	EM372	PKSTTRFC12161620
Stretto MWD 16-20inWx12-16inD Tilt Recline Powerchair	EM372	PKSTTRFC16201216
Stretto MWD 16-20inWx16-20inD Tilt Recline Powerchair	EM372	PKSTTRFC16201620

## Product URL

<https://alliedmedical.co.nz/products/quantum-edge-3-stretto>

# Aviva RWD Powerchair – Moderate Configuration



**Description:** Aviva Rear Wheel Drive (RWD) Narrow or Wide base. Selections of: Tilt, Tilt and elevate, Tilt and recline seating functions. Non-expandable joystick. Optional Centre-mount or 70 Degree Swing away leg rests.

**Cost bracket:** \$13000.00 - \$14000 (Tilt), \$18000.00 - \$19000 (Tilt Elevate or Tilt Recline)

**Supplier:** Invacare New Zealand

**Intended Use:** Power wheelchair with moderate configuration (Adult). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1).

**Key Features:** Package includes: Aviva RX Narrow or Wide Base, ULM seating system (selections of: Tilt, Tilt and elevate, Tilt and recline), 2 point padded hip belt, REM216 non-expandable Joystick, 73amp batteries, centre-mount Leg rest or 70 Degree swing away leg rests, fixed mount joystick, Dual-Post Flip-Back Armrests (Non-Recline Seats) or Cantilever Armrests (Recline seat). (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg

# Aviva RWD Powerchair – Moderate Configuration



**Description:** Aviva Rear Wheel Drive (RWD) Narrow or Wide base. Selections of: Tilt, Tilt and elevate, Tilt and recline seating functions. Non-expandable joystick. Optional Centre-mount or 70 Degree Swing away leg rests.

**Cost bracket:** \$13000.00 - \$14000 (Tilt), \$18000.00 - \$19000 (Tilt Elevate or Tilt Recline)

**Supplier:** Invacare New Zealand

## Product Items on aSAP

Item Name	SIC	Supplier Code
Aviva RWD Nrrw 14-24inWx14-24inD Tilt Elevate Powerchair	EM391	PKG-AVIVA RX NARR 3
Aviva RWD Nrrw 14-24inWx14-24inD Tilt Powerchair	EM389	PKG-AVIVA RX NARR 1
Aviva RWD Nrrw 14-24inWx14-24inD Tilt Recline Powerchair	EM390	PKG-AVIVA RX NARR 2
Item Name	SIC	Supplier Code
Aviva RWD Wide 14-24inWx14-24inD Tilt Elevate Powerchair	EM394	PKG-AVIVA RX WIDE 3
Aviva RWD Wide 14-24inWx14-24inD Tilt Powerchair	EM392	PKG-AVIVA RX WIDE 1
Aviva RWD Wide 14-24inWx14-24inD Tilt Recline Powerchair	EM393	PKG-AVIVA RX WIDE 2

**Note:** The distance between the drive wheel and front casters can be adjusted; for example, longer distances may improve the wheelchair's stability.

## Product URL

<https://www.invacare.co.nz/move/powered-mobility/powered-wheelchairs/aviva-rx>

# Adult Power Wheelchair with High Configuration



Power Wheelchair Categories	Power Wheelchair with High Configuration
Adult	<p><b><u>MWD options:</u></b></p> <p>Q500M (Medifab) TA IQ MWD (Cubro) TDXSP2 REM400 (Invacare) Quantum Edge 3 (Allied Medical) M3 Corpus (Permobil)</p> <p>Q500 Mini (Medifab) Quantum Edge 3 Stretto (Allied Medical)</p> <p><b><u>RWD options:</u></b></p> <p>Aviva RX REM400 (Invacare)</p> <p><b><u>FWD options:</u></b></p> <p>F3 Corpus (Permobil)</p>

Green = power wheelchairs in one category only

Blue = power wheelchairs in more than one category

MWD = Mid Wheel Drive

RWD = Rear Wheel Drive

FWD = Front Wheel Drive

# Q500M MWD Powerchair – High Configuration



**Description:** Q500M MWD Powerchair. Includes Adjustable Seat, Armrests & Footrests. 13in Drive Wheels. 7in Castors. 60 Ah Batteries. 10km Speed. 40km Range. Selections of: 50 Degree Posterior Tilt. 30cm Seat Elevate. 130 Degree Recline. Swing Away Angle Adjustable Footplates or Centre-mount (CM) Footrests.

**Cost bracket:** \$18000.00 - > \$20000

**Supplier:** Medifab Limited

**Intended Use:** Power wheelchair with high configuration (Adult). Suitable for indoor mobility and outdoor urban terrain. Designed for individuals with high and complex needs, including progressive conditions with the potential for changing needs. Features multiple power seat functions, supports alternative controls, and can support individuals displaying WMPM Level 2 clinical indicators.

**Key Features:** Package includes: R-Net Expandable Electronics. 13in Drive Wheels. 7in Castors. Transit Tie Downs. 4 Pole Motors. 60Ah Gel Batteries. Battery Charger. Sedeo Pro Seat Top. Selections of: Power Seat Tilt. Power Seat Elevate. Power Back Recline. Power Swing Away Elevating Legrests or Power Centre Mount Footplate. Cantilever Flip Back Armrests. RHS Swing Away Controller Mount. CJSM1 Standard Joystick. Standard Lap Belt. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg

# Q500M MWD Powerchair – High Configuration



**Description:** Q500M MWD Powerchair. Includes Adjustable Seat, Armrests & Footrests. 13in Drive Wheels. 7in Castors. 60 Ah Batteries. 10km Speed. 40km Range. Selections of: 50 Degree Posterior Tilt. 30cm Seat Elevate. 130 Degree Recline. Swing Away Angle Adjustable Footplates or Centre-mount (CM) Footrests.

**Cost bracket:** \$18000.00 - > \$20000

**Supplier:** Medifab Limited

## Product Items on aSAP

Item Name	SIC	Supplier Code
Q500M MWD 16-20inWx16-22inD Tilt Recline (PCM) Powerchair	EM349	1652-1901-008
Q500M MWD 16-20inWx16-22inD Tilt Recline (PSA) Powerchair	EM347	1652-1901-008

Item Name	SIC	Supplier Code
Q500M MWD 16-20inWx16-22inD TRE (CM) Powerchair	EM346	1652-1901-010
Q500M MWD 16-20inWx16-22inD TRE (PCM) Powerchair	EM350	1652-1901-012
Q500M MWD 16-20inWx16-22inD TRE (PSA) Powerchair	EM348	1652-1901-009
Q500M MWD 16-20inWx16-22inD TRE (SA) Powerchair	EM342	1652-1901-007

## Product URL

<https://medifab.com/nz/product/quickie-q500m/>



# TA iQ MWD Powerchair – High Configuration



**Description:** TA iQ Mid Wheel Drive (MWD) power wheelchair - Tilt/recline/Elevate (TRE) /Leg Elevate (Centre Mount). Base width: 630mm. Base length (without footplates): 910mm. Seat to floor height: 380mm. Turning radius: 450mm. Speed: 12.5km/h. Battery capacity: 80amps. Selections of: Tilt 45°. Backrest recline with mechanical shear reduction. Seat Elevate 11.8in. Powered centre mount (PCM) elevating and articulating legrest with flip up, angle adjustable footplates (2 individual pieces).

**Cost bracket:** \$18000.00 - \$20000

**Supplier:** Cubro

**Intended Use:** Power wheelchair with high configuration (Adult). Suitable for indoor mobility and outdoor urban terrain. Designed for individuals with high and complex needs, including progressive conditions with the potential for changing needs. Features multiple power seat functions, supports alternative controls, and can support individuals displaying WMPM Level 2 clinical indicators.

**Key Features:** Package includes: Standard controller and electronics (R-net expandable). Swing away controller mount / receiver. Centre mount leg hanger with flip up, angle adjustable footplates (2 individual pieces). Flip back, height adjustable arm supports with standard type, full length arm pads. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery charger. (Refer to each product package for specific components included)

**User Weight Capacity:** 140kg

# TA iQ MWD Powerchair – High Configuration



**Description:** TA iQ Mid Wheel Drive (MWD) power wheelchair - Tilt/recline/Elevate (TRE) /Leg Elevate (Centre Mount). Base width: 630mm. Base length (without footplates): 910mm. Seat to floor height: 380mm. Turning radius: 450mm. Speed: 12.5km/h. Battery capacity: 80amps. Selections of: Tilt 45°. Backrest recline with mechanical shear reduction. Seat Elevate 11.8in. Powered centre mount (PCM) elevating and articulating legrest with flip up, angle adjustable footplates (2 individual pieces).

**Cost bracket:** **\$18000.00 - \$20000**

**Supplier:** Cubro

## Product Items on aSAP

Item Name	SIC	Supplier Code
TA iQ MWD 16-24inWx10-24inD TRE (PCM) Powerchair	EM336	XTAIQMWD-MF
TA iQ MWD 16-24inWx10-24inD TRE Powerchair	EM335	XTAIQMWD-TRE

## Product URL

<https://www.cubro.co.nz/healthcare-products/powered-wheelchairs/mid-wheel-drive/ta-iq-mwd-wheelchair>



# TDXSP2 MWD Powerchair – High Configuration



**Description:** TDXSP2 Narrow or Wide Mid wheel drive (MWD) base, Selections of: Tilt, Recline, Elevate (TRE) seating functions. Expandable joystick. Optional Centre-mount, 70 degree Swing away leg rests, or Individual powered swing away (PSA) leg rests.

**Cost bracket:** > \$20000

**Supplier:** Invacare New Zealand

**Intended Use:** Power wheelchair with high configuration (Adult). Suitable for indoor mobility and outdoor urban terrain. Designed for individuals with high and complex needs, including progressive conditions with the potential for changing needs. Features multiple power seat functions, supports alternative controls, and can support individuals displaying WMPM Level 2 clinical indicators.

**Key Features:** Package includes: TDXSP2 Narrow or Wide Base, ULM seating system – selections of Tilt, Recline, Elevate (TRE), 2-point padded hip belt, REM400 expandable Joystick, 52amp batteries, fixed mount joystick. Cantilever Armrests (Recline seat), optional Centre-mount leg rest, 70 degree swing away leg rests, or power elevating leg rest. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg

# TDXSP2 MWD Powerchair – High Configuration



**Description:** TDXSP2 Narrow or Wide Mid wheel drive (MWD) base, Selections of: Tilt, Recline, Elevate (TRE) seating functions. Expandable joystick. Optional Centre-mount, 70 degree Swing away leg rests, or Individual powered swing away (PSA) leg rests.

**Cost bracket:** > \$20000

**Supplier:** Invacare New Zealand

## Product Items on aSAP

Item Name	SIC	Supplier Code
TDXSP2 MWD Nrrw 14-24Wx14-24D Tilt Recline (PSA) Powerchair	EM414	PKG-TDX NARR 5
TDXSP2 MWD Wide 14-24Wx14-24D Tilt Recline (PSA) Powerchair	EM417	PKG-TDX WIDE 5

Item Name	SIC	Supplier Code
TDXSP2 MWD Nrrw 14-24inWx14-24inD TRE Powerchair	EM413	PKG-TDX NARR 4
TDXSP2 MWD Nrrw 14-24Wx14-24inD TRE (PSA) Powerchair	EM415	PKG-TDX NARR 6
TDXSP2 MWD Wide 14-24inWx14-24inD TRE (PSA) Powerchair	EM418	PKG-TDX WIDE 6
TDXSP2 MWD Wide 14-24inWx14-24inD TRE Powerchair	EM416	PKG-TDX WIDE 4

## Product URL

<https://www.invacare.co.nz/move/powerd-mobility/powerd-wheelchairs/tdx-sp2-linx>

# Edge 3 MWD Powerchair – High Configuration



**Description:** A Mid Wheel Drive (MWD) powerchair, easy to customise the mechanical and electronic set up. Ideal for home and in the community. 14in Drive Wheels. 6in Solid Castor Wheels Front. 6in Solid Castor Wheels Rear. Min Seat to Floor Height 16.5in. Maximum Rate Slope 7.5°. Max Speed 10km/h. Selections of: Power tilt, recline, elevate (TRE) and power elevating centre mount (PCM) leg rest.

**Cost bracket:** > \$20000

**Supplier:** Allied Medical Ltd

**Intended Use:** Power wheelchair with high configuration (Adult). Suitable for indoor mobility and outdoor urban terrain. Designed for individuals with high and complex needs, including progressive conditions with the potential for changing needs. Features multiple power seat functions, supports alternative controls, and can support individuals displaying WMPM Level 2 clinical indicators.

**Key Features:** Package includes: Standard controller and electronic (Q-Logic 3 expandable. Swing away controller mount / receiver. Centre mount leg hanger with flip up, angle adjustable footplates. Flip back,. height adjustable arm supports. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery and charger. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg

# Edge 3 MWD Powerchair – High Configuration



## Product Items on aSAP

Item Name	SIC	Supplier Code
Edge 3 MWD 16-20inWx12-16inD Tilt Recline (PCM) Powerchair	EM380	PKE3TRFPC16201216
Edge 3 MWD 16-20inWx16-20inD Tilt Recline (PCM) Powerchair	EM380	PKE3TRFPC16201620
Edge 3 MWD 16-20inWx20-22inD Tilt Recline (PCM) Powerchair	EM380	PKE3TRFPC16202024
Edge 3 MWD 20-22inWx16-20inD Tilt Recline (PCM) Powerchair	EM380	PKE3TRFPC20241620
Edge 3 MWD 20-22inWx20-22inD Tilt Recline (PCM) Powerchair	EM380	PKE3TRFPC20242024

## Product Items on aSAP

Item Name	SIC	Supplier Code
Edge 3 MWD 14-16inWx16-20inD TRE Powerchair	EM379	PKE3TERFC14161620
Edge 3 MWD 16-20inWx16-20inD TRE Powerchair	EM379	PKE3TERFC16201620
Edge 3 MWD 16-20inWx20-22inD TRE Powerchair	EM379	PKE3TERFC16202024
Edge 3 MWD 20-22inWx16-20inD TRE Powerchair	EM379	PKE3TERFC20241620
Edge 3 MWD 20-22inWx20-22inD TRE Powerchair	EM379	PKE3TERFC20242024

Item Name	SIC	Supplier Code
Edge 3 MWD 14-16inWx16-20inD TRE (PCM) Powerchair	EM381	PKE3TERFPC14161620
Edge 3 MWD 16-20inWx16-20inD TRE (PCM) Powerchair	EM381	PKE3TERFPC16201620
Edge 3 MWD 16-20inWx20-22inD TRE (PCM) Powerchair	EM381	PKE3TERFPC16202024
Edge 3 MWD 20-22inWx16-20inD TRE (PCM) Powerchair	EM381	PKE3TERFPC20241620
Edge 3 MWD 20-22inWx20-22inD TRE (PCM) Powerchair	EM381	PKE3TERFPC20242024

## Product URL

<https://alliedmedical.co.nz/products/quantum-edge-3>

# Permobil M3 Powerchair - High Configuration



**Description:** Permobil M3 Mid-Wheel-Drive Power Wheelchair with selections of Power Tilt, Power Elevate, Power Recline and Power Centre-mount (PCM) Elevating Leg rest. Max Speed 10km/h, 14in Drive Wheel, Floor-to-Seat Height - 17.5-19.5in. Selections of: 50° Power Tilt, 12in Power Elevate, 180° Power Recline, 170° Power Elevating Leg rest.

**Cost bracket:** > \$20000

**Supplier:** Permobil New Zealand Limited

**Intended Use:** Power wheelchair with high configuration (Adult). Suitable for indoor mobility and outdoor urban terrain. Designed for individuals with high and complex needs, including progressive conditions with the potential for changing needs. Features multiple power seat functions, supports alternative controls, and can support individuals displaying WMPM Level 2 clinical indicators.

**Key Features:** Package includes: M3 Power Base with Selections of: 50° Power Tilt, 12in Power Elevate, 180° Power Recline, 170° Power Elevating Leg rest, Seat Pan, Width adjustable 17-23in, Depth adjustable 14-22in. R-net Colour CJSM Expandable Joystick, ICS Seating Switchbox, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Anti-Shear Backrest Assembly, Flip Back Height and Width Adjustable Armrests, Manual or Power Elevating Centre-Mount Leg rest, with Calf Supports, USB Device Charger. 10° Active Reach (anterior tilt), Standard 2-point Positioning Belt, Transfer Handles, LED Light Kit, Leatherette or Fabric Armrests. (Refer to each product package for specific components included)

**User Weight Capacity:** 150kg

# Permobil M3 Powerchair - High Configuration



**Description:** Permobil M3 Mid-Wheel-Drive Power Wheelchair with selections of Power Tilt, Power Elevate, Power Recline and Power Centre-mount (PCM) Elevating Leg rest. Max Speed 10km/h, 14in Drive Wheel, Floor-to-Seat Height - 17.5-19.5in. Selections of: 50° Power Tilt, 12in Power Elevate, 180° Power Recline, 170° Power Elevating Leg rest.

**Cost bracket:** > \$20000

**Supplier:** Permobil New Zealand Limited

## Product Items on aSAP

Item Name	SIC	Supplier Code
Permobil M3 17-23inWx14-22inD Tilt Recline (PCM) Powerchair	EM318	PKGM3TRL

Item Name	SIC	Supplier Code
Permobil M3 17-23inWx14-22inD TRE (PCM) Powerchair	EM319	PKGM3TERL
Permobil M3 17-23inWx14-22inD TRE Powerchair	EM316	PKGM3TER

## Product URL

<https://www.permobil.com/en-nz/products/power-wheelchairs/m3-corpus>



# Q500 Mini MWD Powerchair – High Configuration



**Description:** Q500 Mini MWD Powerchair. 13in Drive Wheels. 7in Castors. 50 Ah Batteries. 10km Speed. 30km Range. Selections of: 50 Degree Posterior Tilt. 30cm Seat Elevate. 130 Degree Recline. Power Centre-mount (PCM) Elevating & Articulating Leg rests with Dual Flip Up Footplates or Removable Manual/ Power Swing Away (PSA) Elevating Leg rests with Angle Adjustable Footplates, swings In/Out. Calf Pads Included.

**Cost bracket:** \$19000 - <> \$20000

**Supplier:** Medifab Limited

**Intended Use:** Power wheelchair with high configuration (Adult). Suitable for indoor mobility and outdoor urban terrain. Designed for individuals with high and complex needs, including progressive conditions with the potential for changing needs. Features multiple power seat functions, supports alternative controls, and can support individuals displaying WMPM Level 2 clinical indicators.

**Key Features:** Package includes: R-Net Expandable Electronics. 13in Drive Wheels. 7in Castors. Transit Tie Downs. 4 Pole Motors. 50Ah Gel Batteries. Battery Charger. Selections of: Power Tilt, Power Elevate, Power Recline and Power Elevating Leg rests. Sedeo Pro Seat Top. Dual Flip Up Centre-mount (CM) Footrests or Manual 70 Degree/ Power Elevating Swing Away Leg rests with Aluminium Angle Adjustable Footplates. Cantilever Flip Back Armrests. RHS Swing Away Controller Mount. CJSM1 Standard Joystick. Standard Lap Belt. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg



# Q500 Mini MWD Powerchair – High Configuration



**Description:** Q500 Mini MWD Powerchair. 13in Drive Wheels. 7in Castors. 50 Ah Batteries. 10km Speed. 30km Range. Selections of: 50 Degree Posterior Tilt. 30cm Seat Elevate. 130 Degree Recline. Power Centre-mount (PCM) Elevating & Articulating Leg rests with Dual Flip Up Footplates or Removable Manual/ Power Swing Away (PSA) Elevating Leg rests with Angle Adjustable Footplates, swings In/Out. Calf Pads Included.

**Cost bracket:** \$19000 - <> \$20000

**Supplier:** Medifab Limited

## Product Items on aSAP

Item Name	SIC	Supplier Code
Q500 Mini MWD 16-20Wx16-22inD Tilt Recline (PCM) Powerchair	EM359	1651-1901-008
Q500 Mini MWD 16-20Wx16-22inD Tilt Recline (PSA) Powerchair	EM361	1651-1901-011

Item Name	SIC	Supplier Code
Q500 Mini MWD 16-20inWx16-22inD TRE (CM) Powerchair	EM357	1651-1901-007
Q500 Mini MWD 16-20inWx16-22inD TRE (PCM) Powerchair	EM360	1651-1901-009
Q500 Mini MWD 16-20inWx16-22inD TRE (PSA) Powerchair	EM362	1651-1901-012
Q500 Mini MWD 16-20inWx16-22inD TRE (SA) Powerchair	EM358	1651-1901-010

## Product URL

<https://medifab.com/nz/product/quickie-q500m-mini/>

# Stretto MWD Powerchair – High Configuration



**Description:** The narrowest, most compact and manoeuvrable power base in the Quantum range. 14in Drive Wheels. 5in Solid Castor Wheels Front. 6in Solid Castor Wheels Rear. Min Seat to Floor Height 17.5in. Maximum Rate Slope 6°. Max Speed 10km/h. Mid-Wheel Drive (MWD). Q-Logic 3 Expandable Electronics. SRS Suspension. High Efficiency And Mileage. 2 Pole Motor. Fender lights. Selections of: Power tilt, recline, elevate (TRE), power elevating centre mount (PCM) leg rest.

**Cost bracket:** > \$20000

**Supplier:** Allied Medical Ltd

**Intended Use:** Power wheelchair with high configuration (Adult). Suitable for indoor mobility and outdoor urban terrain. Designed for individuals with high and complex needs, including progressive conditions with the potential for changing needs. Features multiple power seat functions, supports alternative controls, and can support individuals displaying WMPM Level 2 clinical indicators.

**Key Features:** Package includes: Standard controller and electronic (Q-Logic 3 expandable). Swing away controller mount / receiver. Centre mount leg hanger with flip up, angle adjustable footplates. Flip back, height adjustable arm supports. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery and charger. Selections of: Power tilt, recline, elevate (TRE), power elevating centre mount (PCM) leg rest. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg

# Stretto MWD Powerchair – High Configuration



**Description:** The narrowest, most compact and manoeuvrable power base in the Quantum range. 14in Drive Wheels. 5in Solid Castor Wheels Front. 6in Solid Castor Wheels Rear. Min Seat to Floor Height 17.5in. Maximum Rate Slope 6°. Max Speed 10km/h. Mid-Wheel Drive (MWD). Q-Logic 3 Expandable Electronics. SRS Suspension. High Efficiency And Mileage. 2 Pole Motor. Fender lights. Selections of: Power tilt, recline, elevate (TRE), power elevating centre mount (PCM) leg rest.

**Cost bracket:** > \$20000

**Supplier:** Allied Medical Ltd

## Product Items on aSAP

Item Name	SIC	Supplier Code
Stretto MWD 16-20inWx12-16inD Tilt Recline (PCM) Powerchair	EM374	PKSTTRFPC16201216
Stretto MWD 16-20inWx16-20inD Tilt Recline (PCM) Powerchair	EM374	PKSTTRFPC16201620
Stretto MWD 16-20inWx20-22inD Tilt Recline (PCM) Powerchair	EM374	PKSTTRFPC16202024

Item Name	SIC	Supplier Code
Stretto MWD 12-16inWx12-16inD TRE Powerchair	EM373	PKSTTERFC12161216
Stretto MWD 14-16inWx16-20inD TRE (PCM) Powerchair	EM375	PKSTTERFPC14161620
Stretto MWD 14-16inWx16-20inD TRE Powerchair	EM373	PKSTTERFC14161620
Stretto MWD 16-20inWx16-20inD TRE (PCM) Powerchair	EM375	PKSTTERFPC16201620
Stretto MWD 16-20inWx16-20inD TRE Powerchair	EM373	PKSTTERFC16201620

## Product URL

<https://alliedmedical.co.nz/products/quantum-edge-3-stretto>

# Aviva RWD Powerchair – High Configuration



**Description:** Aviva Rear Wheel Drive (RWD) Narrow or Wide base. Selections of: Tilt, Recline, Elevate (TRE) seating functions. Expandable joystick. Optional Centre-mount, 70 degree Swing away leg rests, or Individual powered swing away (PSA) leg rests.

**Cost bracket:** > \$20000

**Supplier:** Invacare New Zealand

**Intended Use:** Power wheelchair with high configuration (Adult). Suitable for indoor mobility and outdoor urban terrain. Designed for individuals with high and complex needs, including progressive conditions with the potential for changing needs. Features multiple power seat functions, supports alternative controls, and can support individuals displaying WMPM Level 2 clinical indicators.

**Key Features:** Package includes: Aviva RX Narrow or Wide Base, ULM seating system – selections of Tilt, Recline, Elevate (TRE), 2-point padded hip belt, REM400 expandable Joystick, 73amp batteries, fixed mount joystick, Dual-Post Flip-Back Armrests (Non-Recline Seats) or Cantilever Armrests (Recline seat), optional Centre-mount leg rest, 70 degree swing away leg rests, or power elevating leg rest. (Refer to each product package for specific components included)

**User Weight Capacity:** 136kg

# Aviva RWD Powerchair – High Configuration



**Description:** Aviva Rear Wheel Drive (RWD) Narrow or Wide base. Selections of: Tilt, Recline, Elevate (TRE) seating functions. Expandable joystick. Optional Centre-mount, 70 degree Swing away leg rests, or Individual powered swing away (PSA) leg rests.

**Cost bracket:** > \$20000

**Supplier:** Invacare New Zealand

## Product Items on aSAP

Item Name	SIC	Supplier Code
Aviva RWD Nrrw 14-24Wx14-24inD Tilt Recline (PSA) Powerchair	EM396	PKG-AVIVA RX NARR 5
Aviva RWD Wide 14-24Wx14-24inD Tilt Recline (PSA) Powerchair	EM399	PKG-AVIVA RX WIDE 5

Item Name	SIC	Supplier Code
Aviva RWD Nrrw 14-24inWx14-24inD TRE (PSA) Powerchair	EM397	PKG-AVIVA RX NARR 6
Aviva RWD Nrrw 14-24inWx14-24inD TRE Powerchair	EM395	PKG-AVIVA RX NARR 4
Aviva RWD Wide 14-24inWx14-24inD TRE (PSA) Powerchair	EM400	PKG-AVIVA RX WIDE 6
Aviva RWD Wide 14-24inWx14-24inD TRE Powerchair	EM398	PKG-AVIVA RX WIDE 4

**Note:** The distance between the drive wheel and front casters can be adjusted; for example, longer distances may improve the wheelchair's stability.

## Product URL

<https://www.invacare.co.nz/move/powering-mobility/powering-wheelchairs/aviva-rx>

# Permobil F3 Powerchair - High Configuration



**Description:** Permobil F3 Mid-Wheel-Drive Power Wheelchair with selections of Power Tilt, Power Elevate, Power Recline and Power Centre-mount (PCM) Elevating Leg rest. Max Speed 10km/h, 14in Drive Wheel, Floor-to-Seat Height - 17.5-19.5in. Selections of: 50° Power Tilt, 12in Power Elevate, 180° Power Recline, 170° Power Elevating Leg rest.

**Cost bracket:** > \$20000

**Supplier:** Permobil New Zealand Limited

**Intended Use:** Power wheelchair with high configuration (Adult). Suitable for indoor mobility and outdoor urban terrain. Designed for individuals with high and complex needs, including progressive conditions with the potential for changing needs. Features multiple power seat functions, supports alternative controls, and can support individuals displaying WMPM Level 2 clinical indicators.

**Key Features:** Package includes: F3 Power Base with Selections of: 50° Power Tilt, 12in Power Elevate, 180° Power Recline, 170° Power Elevating Leg rest, Seat Pan, Width adjustable 17-23in, Depth adjustable 14-22in. R-net Colour CJSM Expandable Joystick, ICS Seating Switchbox, Slimline Retractable Joystick Mount, ESP Tracking Technology, FlexLink adjustable suspension system, 4-point Transit Tie-down hardware, VoltPro charger, Anti-Shear Backrest Assembly, Flip Back Height and Width Adjustable Armrests, Manual or Power Elevating Centre-Mount Leg rest, with Calf Supports, USB Device Charger. 10° Active Reach (anterior tilt), Standard 2-point Positioning Belt, Transfer Handles, LED Light Kit, Leatherette or Fabric Armrests. (Refer to each product package for specific components included)

**User Weight Capacity:** 150kg



# Permobil F3 Powerchair - High Configuration



**Description:** Permobil F3 Mid-Wheel-Drive Power Wheelchair with selections of Power Tilt, Power Elevate, Power Recline and Power Centre-mount (PCM) Elevating Leg rest. Max Speed 10km/h, 14in Drive Wheel, Floor-to-Seat Height - 17.5-19.5in. Selections of: 50° Power Tilt, 12in Power Elevate, 180° Power Recline, 170° Power Elevating Leg rest.

**Cost bracket:** > \$20000

**Supplier:** Permobil New Zealand Limited

## Product Items on aSAP

Item Name	SIC	Supplier Code
Permobil F3 17-23inWx14-22inD Tilt Elevate Powerchair	EM313	PKGf3TE
Permobil F3 17-23inWx14-22inD Tilt Powerchair	EM311	PKGf3T
Permobil F3 17-23inWx14-22inD Tilt Recline (PCM) Powerchair	EM320	PKGf3TRL
Permobil F3 17-23inWx14-22inD Tilt Recline Powerchair	EM315	PKGf3TR
Permobil F3 17-23inWx14-22inD TRE (PCM) Powerchair	EM321	PKGf3TERL
Permobil F3 17-23inWx14-22inD TRE Powerchair	EM317	PKGf3TER

## Product URL

<https://www.permobil.com/en-nz/products/power-wheelchairs/f3-corpus>



# Accessible List

## Bariatric Power Wheelchair Range



# Bariatric Power Wheelchair Products List Overview



Power Wheelchair Categories	Power Wheelchair with Minimal Configuration	Power Wheelchair with Moderate Configuration	Power Wheelchair with High Configuration
Bariatric	-	<p><u>MWD options: (Up to 204kg SWL)</u>  <b>TDXSP2 HD REM216 (Invacare)</b>  <b>Quantum Q6 Edge HD (Allied Medical)</b></p> <p><u>FWD option: (&gt;204kg, up to 272kg)</u>  <b>Quantum Q1450 (Allied Medical)</b></p>	<p><u>MWD options: (Up to 204kg SWL)</u>  <b>TDXSP2 HD REM400 (Invacare)</b>  <b>Quantum Q6 Edge HD (Allied Medical)</b></p> <p><u>FWD option: (&gt;204kg, up to 272kg)</u>  <b>Quantum Q1450 (Allied Medical)</b></p>

**Green** = power wheelchairs in one category only

**Blue** = power wheelchairs in more than one category

MWD = Mid Wheel Drive

RWD = Rear Wheel Drive

FWD = Front Wheel Drive

# Bariatric Power Wheelchair Range

Accessable has defined bariatric power wheelchairs in two sub-categories with the intention of having solutions covering a wide range of users with bariatric considerations.

- **Bariatric options up to 204kg SWL (safe weight limit)**
  - TDXSP2 HD REM216 (Invacare)
  - Quantum Q6 Edge HD (Allied Medical)
- **Bariatric options >204kg and up to 272kg SWL**
  - Quantum Q1450 (Allied Medical)

# Bariatric Power Wheelchair with Moderate Configuration



Power Wheelchair Categories	Power Wheelchair with Moderate Configuration
Bariatric	<p><u>MWD options: (Up to 204kg SWL)</u> TDXSP2 HD REM216 (Invacare) Quantum Q6 Edge HD (Allied Medical)</p> <p><u>FWD option: (&gt;204kg, up to 272kg)</u> Quantum Q1450 (Allied Medical)</p>

Green = power wheelchairs in one category only

Blue = power wheelchairs in more than one category

MWD = Mid Wheel Drive

RWD = Rear Wheel Drive

FWD = Front Wheel Drive

# TDXSP2 HD MWD Powerchair – Moderate Configuration



**Description:** TDXSP2 Heavy Duty Wide Mid wheel drive (MWD) base. Selections of: Tilt, Tilt and elevate, Tilt and recline seating functions. Non-expandable joystick. Optional Centre-mount or 70 Degree Swing away leg rests.

**Cost bracket:** **\$15000.00 - \$16000 (Tilt), > \$20000 (Tilt Elevate or Tilt Recline)**

**Supplier:** Invacare New Zealand

**Intended Use:** Power wheelchair with moderate configuration (Bariatric). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1).

**Key Features:** Package includes: TDXSP2 Heavy Duty Base, ULM seating system (selections of: Tilt, Tilt and elevate, Tilt and recline), 2-point padded hip belt, REM216 non-expandable Joystick, 73amp batteries, Centre-mount Leg rest or 70 Degree swing away leg rests, fixed mount joystick, Dual-Post Flip-Back Armrests (Non-Recline Seats) or Cantilever Armrests (Recline seat). (Refer to each product package for specific components included)

**User Weight Capacity:** 180kg

# TDXSP2 HD MWD Powerchair – Moderate Configuration



**Description:** TDXSP2 Heavy Duty Wide Mid wheel drive (MWD) base. Selections of: Tilt, Tilt and elevate, Tilt and recline seating functions. Non-expandable joystick. Optional Centre-mount or 70 Degree Swing away leg rests.

**Cost bracket:** **\$15000.00 - \$16000 (Tilt), > \$20000 (Tilt Elevate or Tilt Recline)**

**Supplier:** Invacare New Zealand

## Product Items on aSAP

Item Name	SIC	Supplier Code
TDXSP2 HD MWD 14-24inWx14-24inD Tilt Elevate Powerchair	EM402	PKG-TDX HD 2
TDXSP2 HD MWD 14-24inWx14-24inD Tilt Powerchair	EM401	PKG-TDX HD 1
TDXSP2 HD MWD 14-24inWx14-24inD Tilt Recline Powerchair	EM403	PKG-TDX HD 3

## Product URL

<https://www.invacare.co.nz/move/powered-mobility/powered-wheelchairs/tdx-sp2-linx>

<https://www.wheelfreedom.com/products/invacare-tdx-sp2-hd-powerchair/specs>

# Q6 Edge HD MWD Powerchair – Moderate Configuration



**Description:** A heavy duty powerchair engineered to meet the performance needs of the most active user. 14in Drive Wheels. 6in Solid Castor Wheels Front. 6in Solid Castor Wheels Rear. Min Seat to Floor Height 17.75in. Maximum Rate Slope 7.5°. Max Speed 7.5km/h. Selections of: Power tilt. Power tilt and elevate. Power Tilt and Recline.

**Cost bracket:** \$18000.00 - \$19000 (Tilt), \$15000.00 - > \$20000 (Tilt Elevate), >\$20000 (Tilt Recline)

**Supplier:** Allied Medical Ltd

**Intended Use:** Power wheelchair with moderate configuration (Bariatric). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1).

**Key Features:** Package includes: Standard controller and electronic (Q-Logic 3 expandable. Swing away controller mount/receiver. Centre mount leg hanger with flip-up, angle adjustable footplates. Flip back height adjustable arm supports. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery and charger. Selections of: Power tilt. Power tilt and elevate. Power Tilt and Recline. (Refer to each product package for specific components included)

**User Weight Capacity:** 204kg (TB3 Tilt HD, TB3 Recline), 180kg (TB3 Elevate HD with 10in Lift)



# Q6 Edge HD MWD Powerchair – Moderate Configuration



## Product Items on aSAP

Item Name	SIC	Supplier Code
Q6 Edge HD 20-24inWx16-20inD Tilt Powerchair	EM382	PKEHDTFC20241620
Q6 Edge HD 20-24inWx20-24inD Tilt Powerchair	EM382	PKEHDTFC20242024
Q6 Edge HD 24-28inWx16-20inD Tilt Powerchair	EM382	PKEHDTFC24281620
Q6 Edge HD 24-28inWx20-24inD Tilt Powerchair	EM382	PKEHDTFC24282024
Q6 Edge HD 28-32inWx16-20inD Tilt Powerchair	EM382	PKEHDTFC28321620
Q6 Edge HD 28-32inWx20-24inD Tilt Powerchair	EM382	PKEHDTFC28322024

## Product Items on aSAP

Item Name	SIC	Supplier Code
Q6 Edge HD 14-16inWx16-20inD Tilt Elevate Powerchair	EM428	PKEHDTFC14161620
Q6 Edge HD 16-20inWx16-20inD Tilt Elevate Powerchair	EM428	PKEHDTFC16201620
Q6 Edge HD 16-20inWx20-24inD Tilt Elevate Powerchair	EM428	PKEHDTFC16202024
Q6 Edge HD 20-24inWx16-20inD Tilt Elevate Powerchair	EM428	PKEHDTFC20241620
Q6 Edge HD 20-24inWx20-24inD Tilt Elevate Powerchair	EM428	PKEHDTFC20242024
Item Name	SIC	Supplier Code
Q6 Edge HD 20-24inWx16-20inD Tilt Recline Powerchair	EM383	PKEHDTRFC20241620
Q6 Edge HD 20-24inWx20-24inD Tilt Recline Powerchair	EM383	PKEHDTRFC20242024
Q6 Edge HD 24-28inWx16-20inD Tilt Recline Powerchair	EM383	PKEHDTRFC24281620
Q6 Edge HD 24-28inWx20-24inD Tilt Recline Powerchair	EM383	PKEHDTRFC24282024
Q6 Edge HD 28-32inWx16-20inD Tilt Recline Powerchair	EM383	PKEHDTRFC28321620
Q6 Edge HD 28-32inWx20-24inD Tilt Recline Powerchair	EM383	PKEHDTRFC28322024

## Product URL

<https://alliedmedical.co.nz/products/q6-edge-hd>

# Quantum Q1450 Powerchair – Moderate Configuration



**Description:** A bariatric powered mobility solution with the necessary options and sizing to accommodate heavy duty requirements. 14in Drive Wheels. 9in Solid Castor Wheels. Min Seat to Floor Height 16.5in. Maximum Rate Slope 6°. Max Speed 7km/h. Front-Wheel Drive. Q-Logic 3 Expandable Electronics. Swing Arm Suspension. High Efficiency And Mileage. 4 Pole Motors. Heavy Duty Capabilities. Selections of: Power tilt. Tilt and Recline.

**Cost bracket:** >\$20000 (Tilt, Tilt Recline)

**Supplier:** Allied Medical Ltd

**Intended Use:** Power wheelchair with moderate configuration (Bariatric). Suitable for indoor mobility and outdoor use in urban environments. Designed for individuals with moderate to complex needs, presenting a stable condition. This wheelchair is appropriate for a majority of clients with various conditions and can support individuals displaying WMPM Level 1 or Level 2 clinical indicators (note: tilt only for WMPML1).

**Key Features:** Package includes: Standard controller and electronic (Q-Logic 3 expandable. Swing away controller mount/receiver. Centre mount leg hanger with flip-up, angle adjustable footplates. Flip back height adjustable arm supports. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery and charger. Selections of: Power tilt. Tilt and Recline. (Refer to each product package for specific components included)

**User Weight Capacity:** 272kg (TB2 Tilt HD)

# Quantum Q1450 Powerchair – Moderate Configuration



## Product Items on aSAP

Item Name	SIC	Supplier Code
Quantum Q1450 20-24inWx16-20inD (Non-Tilt) Powerchair	EM365	PKQ1450NFC20241620
Quantum Q1450 20-24inWx20-24inD (Non-Tilt) Powerchair	EM365	PKQ1450NFC20242024
Quantum Q1450 24-28inWx16-20inD (Non-Tilt) Powerchair	EM365	PKQ1450NFC24281620
Quantum Q1450 24-28inWx20-24inD (Non-Tilt) Powerchair	EM365	PKQ1450NFC24282024
Quantum Q1450 28-32inWx16-20inD (Non-Tilt) Powerchair	EM365	PKQ1450NFC28321620
Quantum Q1450 28-32inWx20-24inD (Non-Tilt) Powerchair	EM365	PKQ1450NFC28322024

## Product Items on aSAP

Item Name	SIC	Supplier Code
Quantum Q1450 20-24inWx16-20inD Tilt Powerchair	EM367	PKQ1450TFC20241620
Quantum Q1450 20-24inWx20-24inD TB2 Tilt Powerchair	EM366	PKQ1450T2FC20242024
Quantum Q1450 20-24inWx20-24inD Tilt Powerchair	EM367	PKQ1450TFC20242024
Quantum Q1450 24-28inWx16-20inD Tilt Powerchair	EM367	PKQ1450TFC24281620
Quantum Q1450 24-28inWx20-24inD Tilt Powerchair	EM367	PKQ1450TFC24282024
Quantum Q1450 24-28inWx24-28inD TB2 Tilt Powerchair	EM366	PKQ1450T2FC24282428
Quantum Q1450 28-32inWx16-20inD Tilt Powerchair	EM367	PKQ1450TFC28321620
Quantum Q1450 28-32inWx20-24inD Tilt Powerchair	EM367	PKQ1450TFC28322024

Item Name	SIC	Supplier Code
Quantum Q1450 20-24inWx16-20inD Tilt Recline Powerchair	EM368	PKQ1450TRFC20241620
Quantum Q1450 20-24inWx20-24inD Tilt Recline Powerchair	EM368	PKQ1450TRFC20242024
Quantum Q1450 24-28inWx16-20inD Tilt Recline Powerchair	EM368	PKQ1450TRFC24281620
Quantum Q1450 24-28inWx20-24inD Tilt Recline Powerchair	EM368	PKQ1450TRFC24282024
Quantum Q1450 28-32inWx16-20inD Tilt Recline Powerchair	EM368	PKQ1450TRFC28321620
Quantum Q1450 28-32inWx20-24inD Tilt Recline Powerchair	EM368	PKQ1450TRFC28322024

## Product URL

<https://alliedmedical.co.nz/products/quantum-1450-powerchair>

# Bariatric Power Wheelchair with High Configuration



Power Wheelchair Categories	Power Wheelchair with High Configuration
Bariatric	<p><u>MWD options: (Up to 204kg SWL)</u> <b>TDXSP2 HD REM400 (Invacare)</b> <b>Quantum Q6 Edge HD (Allied Medical)</b></p> <p><u>FWD option: (&gt;204kg, up to 272kg)</u> <b>Quantum Q1450 (Allied Medical)</b></p>

**Green** = power wheelchairs in one category only

**Blue** = power wheelchairs in more than one category

MWD = Mid Wheel Drive

RWD = Rear Wheel Drive

FWD = Front Wheel Drive

# TDXSP2 HD MWD Powerchair – High Configuration



**Description:** TDXSP2 Heavy Duty Wide Mid wheel drive (MWD) base. Selections of: Tilt, Recline, Elevate (TRE) seating functions. Expandable joystick. Optional Centre-mount, 70 degree Swing away leg rests, or Individual powered swing away (PSA) leg rests.

**Cost bracket:** > \$20000

**Supplier:** Invacare New Zealand

**Intended Use:** Power wheelchair with high configuration (Bariatric). Suitable for indoor mobility and outdoor urban terrain. Designed for individuals with high and complex needs, including progressive conditions with the potential for changing needs. Features multiple power seat functions, supports alternative controls, and can support individuals displaying WMPM Level 2 clinical indicators.

**Key Features:** Package includes: TDXSP2 Heavy Duty Base, ULM seating system – selections of Tilt, Recline, Elevate (TRE), 2-point padded hip belt, REM400 expandable Joystick, 73amp batteries, fixed mount joystick, Cantilever Armrests (Recline seat), optional Centre-mount leg rest, 70 degree swing away leg rests, or power elevating leg rest. (Refer to each product package for specific components included)

**User Weight Capacity:** 180kg

# TDXSP2 HD MWD Powerchair – High Configuration



**Description:** TDXSP2 Heavy Duty Wide Mid wheel drive (MWD) base. Selections of: Tilt, Recline, Elevate (TRE) seating functions. Expandable joystick. Optional Centre-mount, 70 degree Swing away leg rests, or Individual powered swing away (PSA) leg rests.

**Cost bracket:** > \$20000

**Supplier:** Invacare New Zealand

## Product Items on aSAP

Item Name	SIC	Supplier Code
TDXSP2 HD MWD 14-24Wx14-24inD Tilt Recline (PSA) Powerchair	EM405	PKG-TDX HD 5

Item Name	SIC	Supplier Code
TDXSP2 HD MWD 14-24inWx14-24inD TRE (PSA) Powerchair	EM406	PKG-TDX HD 6
TDXSP2 HD MWD 14-24inWx14-24inD TRE Powerchair	EM404	PKG-TDX HD 4

## Product URL

<https://www.invacare.co.nz/move/power-mobility/power-whe-chairs/tdx-sp2-linx>

<https://www.wheelfreedom.com/products/invacare-tdx-sp2-hd-powerchair/specs>



# Q6 Edge HD MWD Powerchair – High Configuration



**Description:** A heavy duty powerchair engineered to meet the performance needs of the most active user. 14in Drive Wheels. 6in Solid Castor Wheels Front. 6in Solid Castor Wheels Rear. Min Seat to Floor Height 17.75in. Maximum Rate Slope 7.5°. Max Speed 7.5km/h. Selections of: Power tilt, recline, elevate (TRE) and power elevating centre mount (PCM) leg rest.

**Cost bracket:** \$18000.00 - \$19000 (Tilt), \$15000.00 - > \$20000 (Tilt Elevate), >\$20000 (Tilt Recline)

**Supplier:** Allied Medical Ltd

**Intended Use:** Power wheelchair with high configuration (Bariatric). Suitable for indoor mobility and outdoor urban terrain. Designed for individuals with high and complex needs, including progressive conditions with the potential for changing needs. Features multiple power seat functions, supports alternative controls, and can support individuals displaying WMPM Level 2 clinical indicators.

**Key Features:** Package includes: Standard controller and electronic (Q-Logic 3 expandable. Swing away controller mount/receiver. Centre mount leg hanger with flip-up, angle adjustable footplates. Flip back height adjustable arm supports. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery and charger. Selections of: Power tilt, recline, elevate (TRE) and power elevating centre mount (PCM) leg rest. (Refer to each product package for specific components included)

**User Weight Capacity:** 204kg (TB3 Tilt HD, TB3 Recline, PCM elevating leg rest), 180kg (TB3 Elevate HD with 10in Lift)



# Q6 Edge HD MWD Powerchair – High Configuration



## Product Items on aSAP

Item Name	SIC	Supplier Code
Q6 Edge HD 20-24inWx16-20inD Tilt Recline Powerchair (PCM)	EM385	PKEHDTRFPC20241620
Q6 Edge HD 20-24inWx20-24inD Tilt Recline Powerchair (PCM)	EM385	PKEHDTRFPC20242024
Q6 Edge HD 24-28inWx16-20inD Tilt Recline (PCM) Powerchair	EM386	PKEHDTRFPC24281620
Q6 Edge HD 24-28inWx20-24inD Tilt Recline (PCM) Powerchair	EM386	PKEHDTRFPC24282024
Q6 Edge HD 28-32inWx16-20inD Tilt Recline (PCM) Powerchair	EM386	PKEHDTRFPC28321620
Q6 Edge HD 28-32inWx20-24inD Tilt Recline (PCM) Powerchair	EM386	PKEHDTRFPC28322024

## Product Items on aSAP

Item Name	SIC	Supplier Code
Q6 Edge HD 14-16inWx16-20inD TRE (PCM) Powerchair	EM387	PKEHDTERFPC14161620
Q6 Edge HD 14-16inWx16-20inD TRE Powerchair	EM384	PKEHDTERFC14161620
Q6 Edge HD 16-20inWx16-20inD TRE (PCM) Powerchair	EM387	PKEHDTERFPC16201620
Q6 Edge HD 16-20inWx16-20inD TRE Powerchair	EM384	PKEHDTERFC16201620
Q6 Edge HD 16-20inWx20-24inD TRE (PCM) Powerchair	EM387	PKEHDTERFPC16202024
Q6 Edge HD 16-20inWx20-24inD TRE Powerchair	EM384	PKEHDTERFC16202024
Q6 Edge HD 20-24inWx16-20inD TRE (PCM) Powerchair	EM387	PKEHDTERFPC20241620
Q6 Edge HD 20-24inWx16-20inD TRE Powerchair	EM384	PKEHDTERFC20241620
Q6 Edge HD 20-24inWx20-24inD TRE (PCM) Powerchair	EM387	PKEHDTERFPC20242024
Q6 Edge HD 20-24inWx20-24inD TRE Powerchair	EM384	PKEHDTERFC20242024

## Product URL

<https://alliedmedical.co.nz/products/q6-edge-hd>

# Quantum Q1450 Powerchair – High Configuration



**Description:** A bariatric powered mobility solution with the necessary options and sizing to accommodate heavy duty requirements. 14in Drive Wheels. 9in Solid Castor Wheels. Min Seat to Floor Height 16.5in. Maximum Rate Slope 6°. Max Speed 7km/h. Front-Wheel Drive. Q-Logic 3 Expandable Electronics. Swing Arm Suspension. High Efficiency And Mileage. 4 Pole Motors. Heavy Duty Capabilities. Selections of: Power tilt, recline and power elevating centre mount (PCM) leg rest.

**Cost bracket:** >\$20000

**Supplier:** Allied Medical Ltd

**Intended Use:** Power wheelchair with high configuration (Bariatric). Suitable for indoor mobility and outdoor urban terrain. Designed for individuals with high and complex needs, including progressive conditions with the potential for changing needs. Features multiple power seat functions, supports alternative controls, and can support individuals displaying WMPM Level 2 clinical indicators.

**Key Features:** Package includes: Standard controller and electronic (Q-Logic 3 expandable. Swing away controller mount/receiver. Centre mount leg hanger with flip-up, angle adjustable footplates. Flip back height adjustable arm supports. 2 point hip positioning belt. Angle adjustable back canes. Standard tyres. Battery and charger. Selections of: Power tilt, recline and power elevating centre mount (PCM) leg rest. (Refer to each product package for specific components included)

**User Weight Capacity:** 272kg (TB2 Tilt HD)

# Quantum Q1450 Powerchair – High Configuration



**Description:** A bariatric powered mobility solution with the necessary options and sizing to accommodate heavy duty requirements. 14in Drive Wheels. 9in Solid Castor Wheels. Min Seat to Floor Height 16.5in. Maximum Rate Slope 6°. Max Speed 7km/h. Front-Wheel Drive. Q-Logic 3 Expandable Electronics. Swing Arm Suspension. High Efficiency And Mileage. 4 Pole Motors. Heavy Duty Capabilities. Selections of: Power tilt, recline and power elevating centre mount (PCM) leg rest.

**Cost bracket:** >\$20000

**Supplier:** Allied Medical Ltd

## Product Items on aSAP

Item Name	SIC	Supplier Code
Quantum Q1450 20-24Wx16-20inD Tilt Recline (PCM) Powerchair	EM369	PKQ1450TRFPC20241620
Quantum Q1450 20-24Wx20-24inD Tilt Recline (PCM) Powerchair	EM369	PKQ1450TRFPC20242024
Quantum Q1450 24-28Wx16-20inD Tilt Recline (PCM) Powerchair	EM369	PKQ1450TRFPC24281620
Quantum Q1450 24-28Wx20-24inD Tilt Recline (PCM) Powerchair	EM369	PKQ1450TRFPC24282024
Quantum Q1450 28-32Wx16-20inD Tilt Recline (PCM) Powerchair	EM369	PKQ1450TRFPC28321620
Quantum Q1450 28-32Wx20-24inD Tilt Recline (PCM) Powerchair	EM369	PKQ1450TRFPC28322024

## Product URL

<https://alliedmedical.co.nz/products/quantum-1450-powerchair>

# Acknowledgement

## **Evaluation Team**

- Jaynail Comahig
- Jonathan Kwan
- Paul Terry
- Rai Singh
- Sook-Yin Ng

## **Supporting Team/ Managers**

- Aelee Sefo
- Cathy Oliver
- Jonti Collis
- Ngaire Jones
- Tertius van Vuuren

## **Assessor Consultation**

- Mobility Solutions Senior Clinicians